English

**Important Information for Augmentation Patients About  
Mentor® MemoryGel™, Mentor® MemoryGel™ Xtra, and  
Mentor® MemoryGel™ BOOST™  
Silicone Gel-Filled Breast Implants**

October 2023

This page is intentionally left blank.

# Important Information for Augmentation Patients about Mentor® MemoryGel™, Mentor® MemoryGel™ Xtra and Mentor® MemoryGel™ BOOST™ Silicone Gel-Filled Breast Implants

October 2023

<b>Table of Contents</b>	<b>Page No.</b>
GLOSSARY.....	5
1. CONSIDERATIONS FOR SILICONE GEL-FILLED BREAST IMPLANT AUGMENTATION.....	9
1.1. What Gives the Breast Its Shape? .....	9
1.2. What Is a Silicone Gel-Filled Breast Implant?.....	9
1.3. Are You Eligible for Silicone Gel-Filled Breast Implants?.....	10
1.4. Important Factors You Should Consider When Choosing Silicone Gel-Filled Implants .....	10
2. BREAST IMPLANT COMPLICATIONS .....	12
3. MENTOR ROUND GEL CORE STUDY RESULTS FOR AUGMENTATION AND REVISION-AUGMENTATION .....	23
3.1. Overview of the Mentor Round Gel Core Study .....	24
3.2. What Were the 4-Year Follow-up Rates in the Mentor Round Gel Core Study Augmentation Patients? .....	24
3.3. What Were the Benefits for the Mentor Round Gel Core Study Augmentation Patients?.....	24
3.4. What Were the 4-Year Complication Rates in the Mentor Round Gel Core Study Augmentation Patients?.....	25
3.5. What Were the Main Reasons for Reoperation in the Mentor Round Gel Core Study Augmentation Patients?.....	28
3.6. What Were the Main Reasons for Implant Removal in the Mentor Round Gel Core Study Augmentation Patients? .....	29
3.7. What Were Other Clinical Data Findings in the Mentor Round Gel Core Study Augmentation Patients?.....	30
4. MENTOR CPG CORE STUDY RESULTS FOR AUGMENTATION AND REVISION-AUGMENTATION.....	32
4.1. Overview of the Mentor CPG Core Study.....	32
4.2. What Were the 3-Year Follow-up Rates in the Mentor CPG Core Study Augmentation Patients?.....	32
4.3. What Were the Benefits for the Mentor CPG Core Study Augmentation Patients? .....	33
4.4. What Were the 3-Year Complication Rates in the Mentor CPG Core Study Augmentation Patients? .....	34
4.5. What Were the Main Reasons for the Mentor CPG Core Study Reoperation in Augmentation Patients? .....	36
4.6. What Were the Main Reasons for Implant Removal in the Mentor CPG Core Study Augmentation Patients?.....	38
4.7. What Were Other Clinical Data Findings in the Mentor CPG Core Study Augmentation Patients? .....	39
5. SURGERY CONSIDERATIONS FOR RECEIVING BREAST IMPLANTS .....	40
5.1. Surgical Considerations for Breast Augmentation .....	40
5.1.1. What Are the Alternatives to Breast Augmentation with Silicone Breast Implants?.....	40
5.1.2. Choosing a Surgeon.....	40
5.1.3. Implant Shape and Size .....	41
5.1.4. Surface Texturing.....	41

---

5.1.5. Materials Used in the Manufacture of Breast Implants.....	42
5.1.6. Implant Placement.....	42
5.1.7. Incision Sites .....	43
5.1.8. Additional Procedures at the Time of Breast Augmentation.....	43
5.1.9. Palpability .....	43
5.1.10. Surgical Setting and Anesthesia .....	44
5.1.11. Postoperative Care.....	44
5.2. Other Factors to Consider in Revision-Augmentation Surgery.....	44
6. FOLLOW-UP EXAMINATIONS .....	44
6.1. Breast Self-Examinations .....	44
6.2. Recommended Imaging Schedule for Implant RUPTURE Surveillance .....	45
6.3. Mammography .....	45
7. THE TYPES OF SILICONE GEL-FILLED BREAST IMPLANTS AVAILABLE FROM MENTOR .....	45
8. PRODUCT REPLACEMENT POLICY AND LIMITED WARRANTIES .....	45
9. OTHER SOURCES OF ADDITIONAL INFORMATION .....	48
10. PATIENT CHECKLIST – TO BE COMPLETED PRIOR TO SURGERY .....	49
11. MAMMOGRAPHY CENTER FORM .....	55
12. PRIMARY PHYSICIAN FORM .....	57
REFERENCES.....	59

## GLOSSARY

<b>Areola</b>	The pigmented or darker colored area of skin surrounding the nipple of the breast.
<b>Asymmetry</b>	Lack of proportion of shape, size, and/or position between the two breasts.
<b>Autoimmune disease</b>	A disease in which the body mounts an “attack” response to its own tissues or cell types. Normally, the body’s immune mechanism is able to distinguish clearly between what is a normal substance and what is foreign. In autoimmune diseases, this system becomes defective and mounts an attack against normal parts of the body, causing tissue injury. Certain diseases such as rheumatoid arthritis, lupus, and scleroderma are considered to be autoimmune diseases.
<b>Axillary</b>	Pertaining to the armpit area.
<b>Biocompatible</b>	The condition of being compatible with living tissues or systems without being toxic.
<b>Biopsy</b>	The removal and examination of tissues, cells, or fluid from the body.
<b>Body Esteem Scale (BES)</b>	A questionnaire that which asks about a person’s body image.
<b>Breast augmentation</b>	A surgical procedure to increase breast size. For this document, it refers to placement of a breast implant. The first time a breast implant is placed to increase breast size, it is called primary augmentation. All subsequent times the implant is replaced, it is called revision-augmentation.
<b>Breast implant</b>	An internal artificial device or implant intended to replace the breast.
<b>Breast Implant Associated – Anaplastic Large Cell Lymphoma</b>	Breast Implant Associated – Anaplastic Large Cell Lymphoma (BIA-ALCL) is not breast cancer. It is a type of non-Hodgkin’s lymphoma (cancer of the immune system) that can be found in the scar tissue and fluid near the implant but in some cases, it can spread throughout the body.
<b>Breast mass</b>	A lump or body in the breast.
<b>Breast reconstruction</b>	A surgical procedure to replace breast tissue that has been removed due to cancer or trauma or that has failed to develop properly due to a severe breast abnormality.
<b>Calcification</b>	Process of hardening by calcium salts.
<b>Capsular contracture</b>	A tightening of the tissue capsule surrounding an implant, resulting in firmness or hardening of the breast and in squeezing of the implant if severe. Capsular contracture is classified by Baker Grades. Baker Grades III and IV are the most severe. Baker Grade III often results in the need for additional surgery (reoperation) because of pain and possibly abnormal appearance. Baker Grade IV usually results in the need for additional surgery (reoperation) because of pain and unacceptable appearance. Capsular contracture II may also result in the need for additional surgery. Capsular contracture is a known risk for implant rupture. Below is a description of each Baker Grade. <ul style="list-style-type: none"> <li>• Baker Grade I – Normally soft and natural appearance</li> <li>• Baker Grade II – A little firm, but breast looks normal</li> <li>• Baker Grade III – More firm than normal, and looks abnormal (change in shape)</li> <li>• Baker Grade IV – Hard, obvious distortion, and tenderness with pain</li> </ul>
<b>Capsule</b>	Scar tissue that forms around the breast implant. Sometimes this capsule squeezes the implant, resulting in capsular contracture (below).
<b>Capsulectomy</b>	Surgical removal of the scar tissue capsule around the implant.
<b>Capsulorrhaphy</b>	Surgical stitching of a tear in the scar tissue capsule around the implant.
<b>Capsulotomy (closed)</b>	An attempt to break the scar tissue capsule around the implant by pressing or pushing on the outside of the breast. This method does not require surgery but is a known risk for rupture of the implant (see Section 1.4).
<b>Capsulotomy (open)</b>	Surgical incision into the scar tissue capsule around the implant.

---

<b>Congenital anomaly</b>	An abnormal development in part of the body.
<b>Connective tissue disease/disorder (CTD)</b>	A disease, group of diseases, or conditions affecting connective tissue, such as muscles, ligaments, skin, etc. and/or the immune system. Connective tissue diseases (“CTDs”) that involve the immune system include autoimmune diseases such as rheumatoid arthritis, lupus, and scleroderma.
<b>Contour Profile Gel (CPG)</b>	Mentor’s shaped implant. This a tear-drop shaped implant, manufactured with Mentor’s more cohesive (firmer feeling) gel.
<b>Contraindication</b>	A use that is improper and should not be followed. Failure to follow contraindications identified in the labeling could cause serious harm.
<b>Contralateral</b>	Opposite side.
<b>Core Study (Round Gel)</b>	The primary clinical study of augmentation, reconstruction, and revision (revision-augmentation and revision reconstruction) patients that supported the approval of the premarket approval (PMA) application. Safety and effectiveness data are collected yearly through 10 years, with the follow-up from years 4 through 10 being performed as part of a post-approval Round Gel Core Study.
<b>CPG Core Study</b>	The primary clinical study of augmentation, reconstruction, and revision (revision-augmentation and revision reconstruction) patients that supported the approval of the premarket approval (PMA) application. Safety and effectiveness data are collected yearly through 10 years, with the follow-up from years 3 through 10 being performed as part of a post-approval CPG Core Study.
<b>Delayed wound healing</b>	Delayed progress in the healing of an opened wound.
<b>Displacement</b>	Movement of the implant from the usual or proper place.
<b>Epidemiological</b>	Relating to the science of explaining the relationships of factors that determine disease frequency and distribution.
<b>Extracapsular rupture</b>	A type of rupture in which the silicone gel is outside of the scar tissue capsule surrounding the implant.
<b>Extrusion</b>	Skin breakdown with the pressing out of the implant through the surgical wound or skin.
<b>Fibromyalgia</b>	A disorder characterized by chronic pain in the muscles and soft tissues surrounding joints, with tenderness at specific sites in the body. It is often accompanied by fatigue.
<b>Fibrous tissues</b>	Connective tissues composed mostly of fibers.
<b>Granuloma</b>	A lump or mass made of inflammatory cells surrounding a foreign substance due to longstanding inflammation.
<b>Hematoma</b>	A collection of blood within a space.
<b>Hypertrophic scarring</b>	An enlarged scar remaining after the healing of a wound.
<b>Immune response</b>	A bodily response to the presence of a foreign substance.
<b>Infection</b>	Invasion with microorganisms (for example, bacteria, viruses). An infection usually results in fever, swelling, redness, and/or pain.
<b>Inflammation</b>	The response of the body to infection or injury that is characterized by redness, swelling, warmth, pain, and/or loss of function.
<b>Inframammary</b>	Below the breast.
<b>Inframammary fold</b>	The crease at the base of the breast and the chest wall.
<b>Inframammary incision</b>	An incision made in the fold below the breast.
<b>Inpatient surgery</b>	A surgical procedure in which the patient is required to stay overnight in the hospital.

<b>Intracapsular rupture</b>	A type of rupture in which the silicone gel remains inside the scar tissue capsule surrounding the implant.
<b>Lactation</b>	The production and secretion of milk by the breast glands.
<b>Low molecular weight (LMW) silicones</b>	Components of silicone of smaller molecular weight that may bleed (leak) out of silicone gel.
<b>Lymphadenopathy</b>	Enlargement of the lymph node(s).
<b>Malposition</b>	Implant malposition or displacement is when the implant is not in the correct spot in the breast. This could have been due to incorrect placement of the implant during the surgery or due to shifting of the implant position over time.
<b>MRI</b>	Magnetic resonance imaging. A radiographic examination that the ability to detect rupture of silicone gel-filled breast implants.
<b>Mammary</b>	Pertaining to the breast.
<b>Mammography</b>	A type of X-ray examination of the breasts used for detection of cancer.
<b>Mammoplasty</b>	Plastic surgery of the breast.
<b>Mastopexy</b>	Plastic surgery to move sagging breasts into a more elevated position.
<b>Metastatic Disease</b>	Spreading of cancer cells from the original site to other parts of the body.
<b>Migration</b>	Movement of silicone materials outside the breast implant.
<b>Necrosis</b>	Death of cells or tissues.
<b>Outpatient surgery</b>	A surgical procedure in which the patient is not required to stay in the hospital overnight.
<b>Palpate</b>	To feel with the hand.
<b>Palpability</b>	The ability to feel the implant.
<b>Pectoralis</b>	Major muscle of the chest.
<b>Periareolar</b>	Around the darkened or pigmented area surrounding the nipple of the breast.
<b>Plastic surgery</b>	Surgery intended for the improvement of appearance of the body.
<b>Postoperatively</b>	After surgery.
<b>Primary breast augmentation</b>	The first time a breast implant is placed for the purpose of breast augmentation.
<b>Ptosis</b>	Breast sagging that is usually the result of normal aging, pregnancy, or weight loss.
<b>Reoperation</b>	An additional surgery after your first breast implantation.
<b>Revision-Augmentation</b>	Refers to the correction or improvement of a primary augmentation. In the context of this document, it refers to surgical removal and replacement of breast implants that were placed originally for primary breast augmentation.
<b>Rheumatological Disease/ Disorder</b>	A variety of diseases involving connective tissue structures of the body, especially the joints and fibrous tissue. These diseases are often associated with pain, inflammation, stiffness, and/or limitation of motion of the affected parts. Can include autoimmune diseases. Fibromyalgia is a rheumatological disorder.
<b>Rosenberg Self-Esteem Scale</b>	A questionnaire that measures self-esteem.
<b>Rupture</b>	A tear or hole in the implant shell. Silicone implant ruptures may be silent or symptomatic. Ruptures can be intracapsular or extracapsular.
<b>Saline</b>	A solution that is made up of water and a small amount of salt.
<b>Scar revision</b>	A surgical procedure to improve the appearance of a scar.

<b>Seroma</b>	A build-up of the watery portion of the blood in a tissue location.
<b>SF-36 Scale</b>	A questionnaire intended to measure health-related quality of life. It includes questions that measure physical, mental, and social health.
<b>Silicone elastomer</b>	A type of silicone that has elastic properties similar to rubber.
<b>Subglandular placement</b>	Placement of a breast implant underneath and within the breast glands but on top of the chest muscle.
<b>Submuscular placement</b>	Placement of a breast implant wholly or partially underneath the chest muscle.
<b>Surgical incision</b>	A cut made to body tissue during surgery.
<b>Symptom</b>	Any perceptible change in the body or its functions that indicates disease or a phase of a disease.
<b>Symptomatic</b>	Any evidence or sign of disease or disorder reported by the patient.
<b>Symptomatic rupture</b>	A breast implant rupture that is associated with symptoms (such as lumps, persistent pain, swelling, hardening, or change in implant shape). Some silicone breast implant ruptures are symptomatic, but most are silent.
<b>Systemic</b>	Pertaining to or affecting the body as a whole.
<b>Tennessee Self-Concept Scale</b>	A questionnaire that evaluates how the patient sees herself and what she does, likes, and feels.

# IMPORTANT INFORMATION FOR AUGMENTATION PATIENTS ABOUT MENTOR® MEMORYGEL™, MENTOR® MEMORYGEL™ XTRA AND MENTOR® MEMORYGEL™ BOOST™ SILICONE GEL-FILLED IMPLANTS

## 1. Considerations for Silicone Gel-Filled Breast Implant Augmentation

The purpose of this brochure is to help you in making an informed decision about having breast implants for augmentation (breast enlargement) or breast revision-augmentation (replacement) surgery. This brochure is not intended to replace consultation with your surgeon. This educational brochure is set up to provide you with information about risks and benefits of Mentor silicone gel-filled (MemoryGel™) breast implants.


Please read this entire brochure carefully, and if you have any questions or there are things you do not understand, please discuss them with your surgeon before making any decisions. There is a Patient Checklist at the end of this document that summarizes the risks and benefits associated with breast implants. You can use this Checklist as a tool to be sure that you have read and understood the important information in this brochure. Both you and your surgeon can sign the Checklist to document that you understand the risks associated with breast implants and alternative options available to you.

You should wait at least 1-2 weeks after reviewing and considering this information before deciding whether to have primary breast augmentation surgery. In the case of a revision-augmentation; however, your surgeon may find it medically necessary to perform surgery sooner.

### 1.1. What Gives the Breast Its Shape?

The breast consists of milk ducts and glands, surrounded by fatty tissue that provides its shape and feel. The chest muscle (pectoralis major muscle) is located beneath the breast. Factors such as pregnancy (when milk glands are temporarily enlarged), rapid weight loss, and the effects of gravity as you age, combine to stretch the skin, which may cause the breast to droop or sag.

It is important to realize that implants are used to make the breast larger. The implants alone may not adequately lift the breast, or correct the effects of pregnancy, weight loss, or skin stretching. Your surgeon may suggest additional procedures at the time of the breast augmentation, such as mastopexy, to help achieve improved breast lift.

### 1.2. What Is a Silicone Gel-Filled Breast Implant?

A breast implant is a sac (implant shell) of silicone elastomer filled with silicone gel, which is surgically implanted either under your breast tissue or under your chest muscle.

### 1.3. Are You Eligible for Silicone Gel-Filled Breast Implants?

Mentor MemoryGel Silicone Gel-Filled Breast Implants, Mentor MemoryGel Xtra and Mentor MemoryGel BOOST Silicone Gel-Filled Breast Implants are indicated for females for the following uses (procedures):

- **Breast augmentation for women at least 22 years old.** Breast augmentation includes primary breast augmentation to increase the breast size, as well as revision surgery to correct or improve the result of a primary breast augmentation surgery.
- **Breast reconstruction.** Breast reconstruction includes primary reconstruction to replace breast tissue that has been removed due to cancer or trauma or that has failed to develop properly due to a severe breast abnormality. Breast reconstruction also includes revision surgery to correct or improve the result of a primary breast reconstruction surgery. (A separate patient brochure is available for and should be read for breast reconstruction.)

#### Contraindications

Breast implant surgery should not be performed in:

- Women with active infection anywhere in their body.
- Women with existing cancer or pre-cancer of their breast who have not received adequate treatment for those conditions.
- Women who are currently pregnant or nursing.

#### Precautions

Safety and effectiveness have not been established in patients with the following:

- Autoimmune diseases (for example, lupus and scleroderma).
- A weakened immune system (for example, currently taking drugs that weaken the body's natural resistance to disease).
- Conditions that interfere with wound healing and blood clotting.
- Reduced blood supply to breast tissue.
- Radiation to the breast following implantation.
- Clinical diagnosis of depression or other mental health disorders, including body dysmorphic disorder and eating disorders. Please discuss any history of mental health disorders with your surgeon prior to surgery. Patients with a diagnosis of depression, or other mental health disorders, should wait until resolution or stabilization of these conditions prior to undergoing breast implantation surgery.

### 1.4. Important Factors You Should Consider When Choosing Silicone Gel-Filled Implants.

- Breast implants are not lifetime devices, and breast implantation is likely not a one-time surgery. You will likely need additional unplanned surgeries on your breasts because of complications or unacceptable cosmetic outcomes. These additional surgeries can include implant removal with or without replacement, or they can include other surgical procedures. When you have your implants replaced (revision-augmentation), your risk of future complications increases compared to first time (primary) augmentation surgery, so you

should also review the complication rates for revision-augmentation patients to see what future risk rates you may experience.

- Many of the changes to your breast following implantation are irreversible (cannot be undone). If you later choose to have your implant(s) removed and not replaced, you may experience unacceptable dimpling, puckering, wrinkling, or other cosmetic changes of the breast, which can be permanent.
- Breast implants may affect your ability to breastfeed, either by reducing or eliminating milk production.
- Rupture of a silicone gel-filled breast implant is most often silent. This means that neither you nor your surgeon will know that your implants have a rupture most of the time.
- The health consequences of ruptured silicone gel-filled breast implants have not been fully established. Refer to Section 2 (Rupture) below for a review of the literature information regarding this issue.
- If implant rupture is noted by MRI, you should speak to your surgeon about whether you should have the implant removed, with or without replacement.
- With breast implants, routine screening mammography for breast cancer will be more difficult. If you are of the proper age for mammography screening, you should continue to undergo routine mammography screening as recommended by your primary care physician. The implant may interfere with finding breast cancer during mammography. Because the breast and implant are squeezed during mammography, an implant could potentially rupture during the procedure, although this occurrence is rare. More x-ray views are necessary for women with breast implants; therefore, you will receive more exposure to radiation. However, the benefit of having the mammogram to find cancer outweighs the risk of the additional x-rays. Be sure to inform the mammography technologist that you have implants.
- You should perform an examination of your breasts every month for cancer screening; however, this may be more difficult with implants. You should ask your surgeon to help you distinguish the implant from your breast tissue.
- You should perform an examination of your breasts for the presence of lumps, persistent pain, swelling, hardening, or change in implant shape, which may be signs of symptomatic rupture of the implant. These should be reported to your surgeon and possibly evaluated with an ultrasound, mammogram, or MRI to screen for rupture.
- You should inform any other doctor who treats you of the presence of your implants to minimize the risk of damage to the implants.
- Capsular contracture is not to be treated by closed capsulotomy or forceful external compression, which will likely result in implant damage, rupture, folds, and/or hematoma.
- Mentor will continue its ongoing Round Gel and CPG Core Studies through 10 years to further evaluate the long-term safety and effectiveness of these products. In addition, Mentor has initiated a separate, 10-year post-approval study in the U.S. and Canada to address specific issues for which the Round Gel Core Study was not designed to fully answer, as well as to provide a real-world assessment of some endpoints. The endpoints in the large post-approval study include long-term local complications, connective tissue disease (CTD), CTD signs and symptoms, neurological disease, neurological signs and symptoms, offspring issues, reproductive issues, lactation issues, cancer, suicide, mammography issues, and MRI compliance and

results. Mentor will update its labeling as appropriate with the results of the Core and post-approval studies. You should also ask your surgeon if he/she has any available updated clinical information.

- It is important that you read this entire brochure because you need to understand the risks and benefits and to have realistic expectations of the outcome of your surgery.

## **2. Breast Implant Complications**

Undergoing any type of surgical procedure involves risks (some serious) such as the effects of anesthesia, infection, swelling, redness, bleeding, pain, and even death, which need to be balanced against the benefits of the breast augmentation surgery. There are potential complications specific to breast implant surgery and breast implants, as described below. It should also be noted that the cited references include data from augmentation and/or reconstruction patients, as well as from a variety of manufacturers and implant models.

### **WARNING:**

- Breast implants are not considered lifetime devices. The longer people have them, the greater the chances are that they will develop complications, some of which will require more surgery.
- Breast implants have been associated with the development of a cancer of the immune system called Breast Implant-Associated Anaplastic Large Cell Lymphoma (BIA-ALCL). This cancer occurs more commonly in patients with textured breast implants than smooth implants, although rates are not well defined. Some patients have died from BIA-ALCL. Please review detailed information on BIA-ALCL in the Cancer section of this document.
- Patients receiving breast implants have reported a variety of systemic symptoms such as joint pain, muscle aches, confusion, chronic fatigue, and other rheumatology signs and symptoms. Individual patient risk for developing these symptoms has not been well established. Some patients report complete resolution of symptoms when the implants are removed without replacement.

### **• Rupture**

Breast implants are not lifetime devices. Breast implants rupture when the shell develops a tear or hole. Rupture can occur at any time after implantation, but it is more likely to occur the longer the implant is implanted. The following things may cause your implant to rupture: damage by surgical instruments; stressing the implant during implantation and weakening it; folding or wrinkling of the implant shell; excessive force to the chest (for example, during closed capsulotomy; see Section 1.4); trauma; compression during mammographic imaging; and severe capsular contracture. Breast implants may also simply wear out over time. Mentor has done studies to better understand what causes breast implants to rupture or leak gel. These studies might not have identified all causes of rupture and these laboratory studies are ongoing

Silicone gel-filled implant ruptures are most often silent. This means that most of the time neither you nor your plastic surgeon will know if the implant has a tear or hole in the shell. However, sometimes there are symptoms associated with gel implant rupture. These symptoms include hard knots or lumps surrounding the implant or in the armpit, change or loss of size or shape of the breast or implant, pain, tingling, swelling, numbness, burning, or hardening of the breast. If you are experiencing these symptoms, you should discuss them with your physician.

When ultrasound, mammogram, or MRI findings of rupture are found, or if your surgeon determines you have signs or symptoms of rupture, your surgeon may recommend that you have the implant and any gel removed, with or without replacement of the implant. It also may be necessary to remove the tissue capsule as well as the implant, which will involve additional surgery, with associated costs. If you have symptoms such as breast hardness, a change in breast shape or size, and/or breast pain, you should have an ultrasound, mammogram, or MRI to determine whether rupture is present.<sup>1,2</sup>

There are also consequences of rupture. If rupture occurs, silicone gel may either remain within the scar tissue capsule surrounding the implant (intracapsular rupture), move outside the capsule (extracapsular rupture), or gel may move beyond the breast (migrated gel). There is also a possibility that rupture may progress from intracapsular to extracapsular and beyond. There have also been health consequences reported in the literature. See below for details.

### Rupture – Round Gel Core Study

In Mentor's Round Gel Core Study, rupture was assessed for patients who had scheduled MRIs to screen for silent rupture (i.e., part of the MRI cohort) and those who were not assessed for rupture by MRI (i.e., part of the non-MRI cohort). For primary augmentation patients in the MRI cohort, the rupture rate was 2% through 4 years. This means that through 4 years, 2 of every 100 primary augmentation women had at least one ruptured breast implant. There were two primary augmentation patients in the Round Gel Core Study with suspected implant ruptures detected via MRI that were later determined not to be ruptured, and one patient who had suffered a hard fall who was found to have a confirmed rupture when the implant was removed.

For revision-augmentation patients in the MRI cohort, the rupture rate was 7% through 4 years. This means that about 7 of 100 women had at least one ruptured breast implant through 4 years. Three patients had a suspected ruptured implant that was determined not to be ruptured. One woman had removal of her breast implants after MRI, and both implants were ruptured. The other implant ruptures have not yet been confirmed with removal and examination of the implant.

There were no ruptures reported in the non-MRI cohorts for either the primary augmentation or revision-augmentation patients through 4 years. The rupture rate beyond 4 years in Mentor's Round Gel Core Study is currently unknown. The Round Gel Core Study will continue to monitor patients for rupture.

Further information on the estimated incidence rate of rupture for round MemoryGel implants is provided by a limited set of long-term follow-up data from an MRI study in the U.K. (Sharpe and Collis). In this study, textured MemoryGel implants placed subglandularly in 101 patients by a single physician, with follow-up of 4-12 years, were evaluated for rupture status by MRI with confirmation by explantation. Based on their results, at 12 years, the estimated cumulative rate of silent ruptures is 15% for the patients and 9% for implants. By implant, at 12 years, the cumulative rate of 9% is the best statistical estimate, and 19% is a worst-case statistical estimate. By patient, at 12 years, the cumulative rate of 15.1% is the best statistical estimate, and 24.5% is the worst-case statistical estimate. These data are consistent with a published MRI-based rupture study of current silicone gel-filled breast implants from a variety of manufacturers.<sup>3</sup>

## Rupture Rate – CPG Core Study

In Mentor's CPG Core Study, rupture was assessed for patients who had scheduled MRIs to screen for silent rupture (i.e., part of the MRI cohort) and those who were not assessed for rupture by MRI (i.e., part of the non-MRI cohort). For primary augmentation patients in the MRI cohort, the rupture rate was 0% through 3 years. For revision-augmentation patients, the rupture rate was 1% after 3 years. There was one revision-augmentation patient in the Mentor CPG Core Study who was not in the MRI cohort with a ruptured implant that was confirmed with examination of the implant following removal. The rupture rate beyond 3 years in Mentor's CPG Core Study continues to be investigated. The Mentor CPG Core Study will continue to monitor patients for rupture.

## Consequences of Rupture Reported in Literature

Studies of Danish women evaluated with MRI involving a variety of manufacturers and implant models showed that about three-fourths of implant ruptures are intracapsular and the remaining one-fourth are extracapsular.<sup>4</sup> Extracapsular ruptures appear largely to be the result of closed capsulotomy (see Section 1.4) and/or trauma to the chest area. For example, there was a significantly higher prevalence of extracapsular ruptures (14.7%) in these Danish women who had undergone closed capsulotomy as compared to those who had not.<sup>5</sup> In a study of British women, one patient observed to have severe bilateral silicone granulomas and bilateral extracapsular ruptures suffered a fractured sternum in a traffic accident.<sup>6</sup>

There is a possibility that rupture may progress from intracapsular to extracapsular and beyond. Studies of Danish women indicate that over a 2-year period, about 10% of the implants with intracapsular rupture progressed to extracapsular rupture as detected by MRI.<sup>7</sup> This means that for women with silicone gel rupture within the scar tissue capsule detected via MRI after 2 years, 1 in 10 of these women will have progression of the gel outside the scar tissue capsule. Approximately half of the women whose ruptures had progressed from intra- to extracapsular reported that they experienced trauma to the affected breast during this time period or had undergone mammography. In the other half, no cause was given. In the women with extracapsular rupture, after 2 years, the amount of silicone seepage outside the scar tissue capsule increased for about 14% of these women. This means that for 100 women with silicone gel rupture outside the scar tissue capsule, the amount of gel outside the scar tissue capsule increased for 14 women 2 years later. This type of information pertains to a variety of silicone implants from a variety of manufacturers and implant models, and it is not specific to Mentor implants.

Below is a summary of information related to the health consequences of implant rupture, which have not been fully established. These reports were in women who had implants from a variety of manufacturers and implant models.

- Local breast complications reported in the published literature that were associated with rupture include breast hardness, a change in breast shape or size, and breast pain.<sup>8</sup> These symptoms are not specific to rupture, as they also are experienced by women who have capsular contracture. Most of the Danish women evaluated in these studies, whose ruptured implants were left in place for two years, reported no symptoms.

- There have been rare reports of gel movement to nearby tissues such as the chest wall, armpit, or upper abdominal wall, and to more distant locations down the arm or into the groin. This has led to nerve damage, granuloma formation (see glossary) and/or breakdown of tissues in direct contact with the gel in a few cases. There have been reports of silicone presence in the liver of patients with silicone breast implants. Movement of silicone gel material to lymph nodes in the axilla also has been reported, even in women without evidence of rupture, leading to lymphadenopathy.<sup>9</sup>
- Concerns have been raised over whether ruptured implants are associated with the development of connective tissue or rheumatic diseases and/or symptoms such as fatigue and fibromyalgia.<sup>10,11,12,13</sup> A number of epidemiology studies have evaluated large populations of women with breast implants from a variety of manufacturers and implant models. These studies do not, taken together, support a significant association of breast implants with a typical, diagnosed rheumatic disease.<sup>14</sup> Other than one small study,<sup>15</sup> few studies distinguish whether the women had ruptured or intact implants.

- **Capsular Contracture**

The scar tissue (capsule) that normally forms around the implant may tighten over time and compress the implant, making it feel firm and leading to what is called capsular contracture. Capsular contracture may be more common following infection, hematoma, and seroma, and the chance of it happening may increase over time. Capsular contracture occurs more commonly in revision-augmentation than in primary augmentation. Because you may have your initial implants replaced, you should be aware that your risk of capsular contracture increases with revision-augmentation. Capsular contracture is a risk factor for implant rupture, and it is one of the most common reasons for reoperation.

Symptoms of capsular contracture range from mild firmness and mild discomfort to severe pain, distorted shape of the implant, and palpability (ability to feel the implant). Capsular contracture is graded into 4 levels depending on its severity:

Baker Grade I:	the breast is normally soft and looks natural
Baker Grade II:	the breast is a little firm but looks normal
Baker Grade III:	the breast is firm and looks abnormal
Baker Grade IV:	the breast is hard, painful, and looks abnormal

Baker Grades III and IV are considered severe and often additional surgery is needed to correct these grades.

Additional surgery may be needed in cases where pain and/or firmness are severe. This surgery ranges from removal of the implant capsule tissue, to removal and possible replacement of the implant itself. This surgery may result in loss of your breast tissue. Capsular contracture may happen again after these additional surgeries. Capsular contracture may increase the risk of rupture.<sup>16</sup>

#### Capsular Contracture – Round Gel Core Study

In Mentor's Round Gel Core Study, for women receiving augmentation implants for the first time, the risk of severe capsular contracture was 9% through 4 years. This means that 9 out of every 100 women who

received Mentor implants for primary breast augmentation had severe capsular contracture at least once during the first 4 years after receiving the implants.

For women receiving revision-augmentation implants, the risk of severe capsular contracture was 20% through 4 years. This means that 20 out of every 100 women who received Mentor implants for breast revision-augmentation had severe capsular contracture at least once during the first 4 years after receiving the implants.

#### Capsular Contracture – CPG Core Study

In Mentor's CPG Core Study, for women receiving augmentation implants for the first time, the risk of severe capsular contracture was 1% through 3 years. This means that approximately 1 out of every 100 women who received Mentor CPG breast implants for primary breast augmentation had severe capsular contracture at least once during the first 3 years after receiving the implants.

For women receiving revision-augmentation implants, the risk of severe capsular contracture was 4% through 3 years. This means that approximately 4 out of every 100 women who received Mentor CPG breast implants for breast revision-augmentation had severe capsular contracture at least once during the first 3 years after receiving the implants.

#### • **Additional Surgeries (Reoperations) – Round Gel Core Study**

You should assume that you will need to have additional surgeries (reoperations). In the Mentor Round Gel Core Study, the reoperation rate was 16% for primary augmentation patients, which means that 16 out of every 100 women who received Mentor implants for primary augmentation had a reoperation during the first 4 years after receiving the implants. The reoperation rate was 28% for revision-augmentation patients, which means that approximately 28 out of every 100 women who received Mentor implants for revision-augmentation had a reoperation during the first 4 years after receiving the implants.

Patients may decide to change the size or type of their implants, requiring additional surgery. Problems such as rupture, capsular contracture, hypertrophic scarring (irregular, raised scar), asymmetry, infection, and shifting can require additional surgery. Summary tables are provided in Section 3.5 that describe the reasons for performing additional surgeries in the Round Gel Core Study. For women receiving primary augmentation implants, the three most common reasons for reoperation were severe capsular contracture, followed by patient request for size/style change, followed by scarring/hypertrophic scarring and hematoma/seroma. For women receiving revision-augmentation implants, the three most common reasons for additional surgery were severe capsular contracture, followed by patient request for style/size change, followed by biopsy.

#### • **Additional Surgeries (Reoperations) – CPG Core Study**

You should assume that you will need to have additional surgeries (reoperations). In the Mentor CPG Core Study, the reoperation rate was 13% for primary augmentation patients, which means that about 13 out of every 100 women who received Mentor CPG breast implants for primary augmentation had a reoperation during the first 3 years after receiving the implants. The reoperation rate was 15% for revision-augmentation

patients, which means that about 15 out of every 100 women who received Mentor CPG breast implants for revision-augmentation had a reoperation during the first 3 years after receiving the implants.

Patients may decide to change the size or type of their implants, requiring additional surgery. Problems such as rupture, capsular contracture, hypertrophic scarring (irregular, raised scar), asymmetry, infection, and shifting can require additional surgery. Summary tables are provided in Section 3.5 that describe the reasons for performing additional surgeries in the CPG Core Study. For women receiving primary augmentation implants, the three most common reasons for reoperation were patient request for size/style change, breast mass, and ptosis. For women receiving revision-augmentation implants, the three most common reasons for additional surgery were wound dehiscence, patient request for size/style change, position dissatisfaction and ptosis.

### • **Implant Removal**

Because these are not lifetime devices, the longer you have your implants, the more likely it will be for you have them removed for any reason, either because of dissatisfaction, an unacceptable cosmetic result, or a complication such as severe capsular contracture. Having your implants removed and replaced increases your chances of getting future complications.

Most women who have their implants removed, have them replaced with new implants, but some women do not. If you choose not to replace your implants, you may have cosmetically unacceptable dimpling, puckering, wrinkling, and/or other potentially permanent cosmetic changes of the breast following removal of the implant. Even if you have your implants replaced, implant removal may result in loss of your breast tissue. Also, implant replacement increases your risks of future complications. For example, the risks of severe capsular contracture and reoperation double for patients with implant replacement compared to first time placement. You should consider the possibility of having your implants replaced and its consequences when making your decision to have implants.

#### Implant Removal - Round Gel Core Study

For women receiving primary augmentation implants in Mentor's Round Gel Core Study, 5% had their implants removed at least once through 4 years. For women receiving revision-augmentation implants in Mentor's Round Gel Core Study, 13% had their implants removed at least once through 4 years. Patient request for style/size change and severe capsular contracture were the most common reasons for implant removal.

#### Implant Removal - CPG Core Study

For women receiving primary augmentation implants in Mentor's CPG Core Study, 5% had their implants removed at least once through 3 years. Patient request for size change and asymmetry were the most common reasons for implant removal. For women receiving revision-augmentation implants in Mentor's CPG Core Study, 11% had their implants removed at least once through 3 years. The most common reason was patient request for size change.

- **Unsatisfactory Results**

Unsatisfactory results such as wrinkling, asymmetry, implant displacement (shifting), incorrect size, unanticipated shape, implant palpability, scar deformity, and/or hypertrophic scarring, may occur. Some of these results may cause discomfort. Pre-existing asymmetry may not be entirely correctable by implant surgery. Revision surgery may be recommended to maintain patient satisfaction but carries additional considerations and risks. Selecting an experienced plastic surgeon may minimize, but not necessarily prevent, unsatisfactory results.

- **Pain**

Pain of varying intensity and length of time may occur and persist following breast implant surgery. In addition, improper size, placement, surgical technique, or capsular contracture may result in pain. You should tell your surgeon about significant pain or if your pain persists.

- **Changes in Nipple and Breast Sensation**

Feeling in the nipple and breast can increase or decrease after implant surgery. The range of changes varies from intense sensitivity to no feeling in the nipple or breast following surgery. While some of these changes can be temporary, they can also be permanent, and may affect your sexual response or your ability to nurse a baby. (See the paragraph on breastfeeding below.)

- **Infection**

Infection can occur with any surgery or implant. Most infections resulting from surgery appear within a few days to weeks after the operation. However, infection is possible at any time after surgery. In addition, breast and nipple piercing procedures may increase the possibility of infection. Infections in tissue with an implant present are harder to treat than infections in tissue without an implant. If an infection does not respond to antibiotics, the implant may have to be removed, and another implant may be placed after the infection is resolved (cleared up). As with many other surgical procedures, in rare instances, toxic shock syndrome has been noted in women after breast implant surgery, and it is a life-threatening condition. Symptoms include sudden fever, vomiting, diarrhea, a sunburn-like rash, red eyes, dizziness, lightheadedness, muscles aches, and/or drops in blood pressure, which may cause fainting. You should contact your doctor immediately for diagnosis and treatment if you have these symptoms.

- **Hematoma/Seroma**

Hematoma is a collection of blood within the space around the implant, and seroma is a build-up of fluid around the implant. Having a hematoma and/or seroma following surgery may result in infection and/or capsular contracture later on. Symptoms from a hematoma or seroma may include swelling, pain, and bruising. If a hematoma or seroma occurs, it will usually be soon after surgery. However, this can also occur at any time after injury to the breast. While the body absorbs small hematomas and seromas, some will require surgery, typically involving draining and potentially placing a surgical drain in the wound temporarily for proper healing. A small scar can result from surgical draining. Implant rupture also can occur from surgical draining if there is damage to the implant during the draining procedure.

- **Breastfeeding**

Breastfeeding difficulties have been reported following breast surgery, including breast reduction and breast augmentation. If your surgeon uses a periareolar surgical approach (an incision around the colored portion surrounding the nipple), it may further increase the chance of breastfeeding difficulties.

- **Calcium Deposits in the Tissue Around the Implant**

Calcium deposits can form in the tissue capsule surrounding the implant. Symptoms may include pain and firmness. Deposits of calcium can be seen on mammograms and can be mistaken for possible cancer, resulting in additional surgery for biopsy and/or removal of the implant to distinguish calcium deposits from cancer. If additional surgery is necessary to examine and/or remove calcifications, this may cause damage to the implants. Calcium deposits also occur in women who undergo breast reduction procedures, in patients who have had hematoma formation, and even in the breasts of women who have not undergone any breast surgery. The occurrence of calcium deposits increases significantly with age.

- **Extrusion**

Extrusion is when the breast implant comes through your skin. This may occur, for example, when your wound has not closed or when breast tissue covering your implants weakens. Radiation therapy has been reported to increase the likelihood of extrusion. Extrusion requires additional surgery and possible removal of the implant, which may result in additional scarring and/or loss of your breast tissue.

- **Necrosis**

Necrosis is the death of cells or tissues. This may prevent or delay wound healing and require surgical correction, which may result in additional scarring and/or loss of your breast tissue. Implant removal may also be necessary. Factors associated with increased necrosis include infection, use of steroids, smoking, chemotherapy, radiation, and excessive heat or cold therapy.

- **Delayed Wound Healing**

Some patients may experience a prolonged wound healing time. Delayed wound healing may increase the risk of infection, extrusion, and necrosis. Depending on the type of surgery or the incision, wound healing times may vary. Smoking may interfere with the healing process. You should contact your surgeon immediately if your wound does not heal within the period of time he/she has discussed with you.

- **Breast Tissue Atrophy/Chest Wall Deformity**

The pressure of the breast implant may cause breast tissue thinning (with increased implant visibility and palpability) and chest wall deformity. This can occur while implants are still in place or following implant removal without replacement. Either of these conditions may result in additional surgeries and/or unacceptable dimpling/puckering of the breast.

---

- **Lymphadenopathy**

Lymphadenopathy is a chronic enlargement of the lymph nodes. A lymph node is a round mass of tissue that makes cells as part of your immune system. The lymph nodes in the armpit (axilla) drain the breast area of fluid. Sometimes the enlarged lymph nodes are painful. If they become too large or painful, the lymph node(s) may need to be surgically removed. Painful and/or enlarged lymph nodes should be reported to your doctor.

Literature reports associate lymphadenopathy with both intact and ruptured silicone breast implants. One study reported that armpit lymph nodes from women with both intact and ruptured silicone gel implants had abnormal tissue reactions, granulomas, and the presence of silicone.<sup>17</sup> These reports were in women who had implants from a variety of manufacturers and implant models.

- **Other Reported Conditions**

There have been reports in the literature of other conditions in women with silicone gel-filled breast implants. Many of these conditions have been studied to evaluate their potential association with breast implants. Although no cause and effect relationship has been established between breast implants and the conditions listed below, you should be aware of these reports. Furthermore, there is the possibility of risks, yet unknown, which in the future could be determined to be associated with breast implants.

- **Connective Tissue Disease (CTD)**

Connective tissue diseases include diseases such as lupus, scleroderma, and rheumatoid arthritis. Fibromyalgia, another CTD, is a disorder characterized by chronic pain in the muscles and soft tissues surrounding joints, with tenderness at specific sites in the body. It is often accompanied by fatigue. There have been a number of published epidemiological studies that have looked at whether having a breast implant is associated with having a typical or defined connective tissue disease. The study size needed to conclusively rule out a smaller risk of connective tissue disease among women with silicone gel-filled breast implants would need to be very large.<sup>18,19,20,21,22,23,24,25,26,27</sup> The published studies taken together show that breast implants are not significantly associated with a risk of developing a typical or defined connective tissue disease.<sup>28,29,30,31</sup> These studies do not distinguish between women with intact and ruptured implants. Only one study evaluated specific connective tissue disease diagnoses and symptoms in women with silent ruptured versus intact implants, but it was too small to rule out a small risk.<sup>32</sup>

- **CTD Signs and Symptoms**

Literature reports have also been made associating silicone breast implants with various rheumatological signs and symptoms such as fatigue, exhaustion, joint pain and swelling, muscle pain and cramping, tingling, numbness, weakness, and skin rashes. Scientific expert panels and literature reports have found no evidence of a consistent pattern of signs and symptoms in women with silicone breast implants.<sup>33,34,35,36,37</sup> Having these rheumatological signs and symptoms does not necessarily mean you have a connective tissue disease; however, you should be aware that you may experience these signs and symptoms after undergoing breast implantation. If you notice an increase in these signs or symptoms, you should consider seeing a

rheumatologist to determine whether these signs or symptoms are due to a connective tissue disorder or autoimmune disease.

- **Immunotoxicity**

While there is no scientific evidence that silicone can cause hypersensitivity reactions in humans, some animal testing reports in the literature suggest that silicone gel may cause an adjuvant effect. The biological mechanism and clinical significance for these findings in animal models remain unknown.

- **Cancer**

Breast Cancer – Reports in the medical literature indicate that patients with breast implants are not at a greater risk than those without breast implants for developing breast cancer.<sup>38,39,40,41,42</sup> Some reports have suggested that breast implants may interfere with or delay breast cancer detection by mammography and/or biopsy; however, other reports in the published medical literature indicate that breast implants neither significantly delay breast cancer detection nor adversely affect cancer survival of women with breast implants.<sup>43,44,45,46,47,48</sup> A review of the available published medical literature concludes that women with breast implants do not have an increased risk of breast cancer or death associated with breast cancer, and that, if they do have breast cancer, the detection of their cancer and their survival are not affected by their implants.<sup>49</sup>

Brain cancer – One study has reported an increased incidence of brain cancer in women with breast implants as compared to the general population.<sup>50</sup> The incidence of brain cancer, however, was not significantly increased in women with breast implants when compared to women who had other plastic surgeries. Other published reviews of large studies in women with cosmetic implants concluded that the evidence does not support an association between brain cancer and breast implants.<sup>51,52</sup>

Respiratory/lung cancer – One study has reported an increased incidence of respiratory/lung cancer in women with breast implants.<sup>53</sup> Other studies of women in Sweden and Denmark have found that women who get breast implants are more likely to be current smokers than women who get breast reduction surgery or other types of cosmetic surgery.<sup>54,55,56</sup>

Cervical/vulvar cancer – One study has reported an increased incidence of cervical/vulvar cancer in women with breast implants.<sup>57</sup> The cause of this increase is unknown, but it has been suggested that the increased risk observed is more likely to reflect other factors in these patients that put them at higher risk of getting this kind of cancer (such as multiple sexual partners).<sup>58</sup>

Lymphomas, including Breast Implant Associated Anaplastic Large Cell Lymphoma (BIA-ALCL)

If you have breast implants, you have a very small, but increased risk of developing breast implant associated anaplastic large cell lymphoma, or BIA-ALCL. BIA-ALCL is not breast cancer — it is a rare type of non-Hodgkin's lymphoma (cancer of the immune system). In most cases, BIA-ALCL is found in the scar tissue and fluid near the implant, but in some cases, it can spread throughout the body. In the cases that have spread beyond the scar tissue and fluid near the implant, rare cases of death have been reported.

Most patients were diagnosed with BIA-ALCL when they sought medical treatment for implant-related symptoms such as swelling, pain, lumps, or asymmetry that developed after their initial surgical sites were

fully healed. In the cases known to FDA to date, BIA-ALCL was diagnosed years after the breast implant was placed. The earliest report was one year after implant placement and the latest was 23 years after the implant surgery. About half the cases occurred within the first 7 years after implant. BIA-ALCL was most often diagnosed in women who had textured implants. The textured implant may have been placed at the most recent surgery or at any other prior breast implant operation. Reports in the medical literature show that high-surface-area textured breast implants are associated with an increased risk of developing BIA-ALCL as compared to low surface-area textured implants.<sup>59,60,61</sup>

If you develop swelling or pain around your breast implants, be sure to talk to your health care provider. Your health care provider should consider the possibility of BIA-ALCL if after you have recovered from your breast implant operation, you later notice changes in the way your breast looks or feels — including swelling or pain around the implant. If your health care provider suspects BIA-ALCL, they will refer you to an appropriate specialist for evaluation which may involve obtaining fluid and tissue samples from around your breast implant. If a diagnosis of BIA-ALCL is confirmed, the doctor will develop an individualized treatment plan for you. Because of the small number of cases worldwide and the variety of available treatment options, there is no single defined treatment. However, if you are diagnosed with BIA-ALCL, the National Comprehensive Cancer Network (NCCN) recommends removing the implant and the surrounding tissue.

If you have breast implants, you should monitor them and follow your routine medical care. You do not need to take any additional steps. It is not necessary to remove your breast implants if you have no symptoms and you have not been diagnosed with BIA-ALCL.

Other cancers – One study has reported an increased incidence of stomach cancer and leukemia in women with breast implants compared to the general population.<sup>62</sup> This increase was not significant when compared to women who had other types of plastic surgeries.

### • **Neurological Disease, Signs, and Symptoms**

Some women with breast implants have complained of neurological symptoms (such as difficulties with vision, sensation, muscle strength, walking, balance, thinking or remembering things) or diseases (such as multiple sclerosis), which they believe are related to their implants. A scientific expert panel report found that the evidence for a neurological disease or syndrome caused by or associated with breast implants is insufficient or flawed.<sup>63</sup>

### • **Suicide**

In several studies, a higher incidence of suicide was observed in women with breast implants.<sup>64,65,66,67</sup> The reason for the observed increase is unknown, but it was found that women with breast implants had higher rates of hospital admission due to psychiatric causes prior to surgery, as compared with women who had breast reduction or in the general population of Danish women.<sup>68</sup>

### • **Effects on Children**

At this time, it is not known if a small amount of silicone may pass through from the breast implant silicone shell into breast milk during breastfeeding. Although there are no current established methods for accurately

detecting silicone levels in breast milk, a study measuring silicon (one component in silicone) levels did not indicate higher levels in breast milk from women with silicone gel-filled implants when compared to women without implants.<sup>69</sup>

In addition, concerns have been raised regarding potential damaging effects on children born to mothers with implants. Two studies in humans have found that the risk of birth defects overall is not increased in children born after breast implant surgery.<sup>70,71</sup> Although low birth weight was reported in a third study, other factors (for example, lower pre-pregnancy weight) may explain this finding.<sup>72</sup> This author recommended further research on infant health. A review of the available literature on the potential effects on children born to mothers with breast implants concludes that the evidence does not support an association between breast implants and health problems in children.<sup>73</sup>

#### • **Potential Health Consequences of Gel Bleed**

Small quantities of low molecular weight (LMW) silicone compounds, as well as platinum (in zero oxidation state), have been found to diffuse (“bleed”) through an intact implant shell.<sup>74,75</sup> Studies on implants implanted for a long duration have suggested that such bleed may be a contributing factor in the development of capsular contracture<sup>76</sup> and lymphadenopathy.<sup>77</sup> Evidence against gel bleed being a significant contributing factor to capsular contracture and other local complications is provided by the fact that there are similar or lower complication rates for silicone gel-filled breast implants than for saline-filled breast implants. Saline-filled breast implants do not contain silicone gel and, therefore, gel bleed is not an issue for those products. Furthermore, toxicology testing has indicated that the silicone material used in the Mentor implants does not cause toxic reactions when large amounts are administered to test animals. It should also be noted that studies reported in the literature have demonstrated that the low concentration of platinum contained in breast implants is in the zero oxidation (most biocompatible) state.<sup>78</sup> In addition, two separate studies sponsored by Mentor have demonstrated that the low concentration of platinum contained in its breast implants is in the zero oxidation (most biocompatible) state.

Mentor performed a laboratory test to analyze the silicones and platinum (used in the manufacturing process), which may bleed out of intact implants into the body. Over 99% of the LMW silicones and platinum stayed in the implant. The overall body of available evidence supports that the extremely low level of gel bleed is of no clinical consequence.

### **3. Mentor Round Gel Core Study Results for Augmentation and Revision-Augmentation**

This section of this brochure summarizes the results of the Mentor Round Gel Core Study conducted on Mentor’s round silicone gel-filled breast implants for primary augmentation and revision-augmentation. The Mentor Round Gel Core Study is the primary clinical study for this product. The results of the Round Gel Core Study give you useful information on the experience of other women with Mentor silicone gel-filled implants. While the results cannot be used to predict your individual outcome, they can be used as a rough guide of what you may expect. Your own complications and benefits depend on many individual factors.

As a note, supplemental safety information was also obtained from the Mentor Adjunct Study, the U.K. Sharpe/Collis Study, and the literature to help assess long-term rupture rate and the consequences of rupture for

this product. The literature, which had the most available information on the consequences of rupture, was also used to assess other potential complications associated with silicone gel-filled breast implants. The key literature information is referenced throughout the Breast Implant Complications section above.

### **3.1. Overview of the Mentor Round Gel Core Study**

The Mentor Round Gel Core Study is a 10-year study to assess safety and effectiveness in augmentation, reconstruction, and revision (revision-augmentation and revision-reconstruction) patients. Patient follow-up is at 6 months and annually through 10 years. Safety is assessed by complications, such as implant rupture, capsular contracture, and reoperation. Benefit (effectiveness) is assessed by patient satisfaction and measures of quality of life (QoL).

The Mentor Round Gel Core Study consists of 1,008 patients, including 552 primary augmentation patients, 145 revision-augmentation patients, 251 primary reconstruction patients, and 60 revision-reconstruction patients. Of these patients, 202 primary augmentation patients, 56 revision-augmentation patients, 134 primary reconstruction patients, and 28 revision-reconstruction patients are in the MRI cohort, which means that they are assessed for silent rupture by MRI at years 1, 2, 4, 6, 8, and 10. The study is currently ongoing, with the results through 4 years reported in this brochure. Mentor will periodically update this brochure as more information becomes available. You should also ask your surgeon if he/she has any available updated clinical information.

Mentor's Round Gel Core Study results indicate that the risk of at least one occurrence of any complication (including reoperation) at some point through 4 years after implant surgery is 39% for primary augmentation patients and 50% for revision-augmentation patients. The information below provides more details about the complications and benefits you may experience.

Described below are the benefits and complications reported in the Mentor Round Gel Core Study for augmentation patients. The findings are described separately for primary augmentation and revision-augmentation patients.

### **3.2. What Were the 4-Year Follow-up Rates in the Mentor Round Gel Core Study Augmentation Patients?**

At the 4-year follow-up visit, data are reported for 80% of the eligible primary augmentation patients and 79% of the eligible revision-augmentation patients.

### **3.3. What Were the Benefits for the Mentor Round Gel Core Study Augmentation Patients?**

The Round Gel Core Study measured a variety of outcomes that assessed the benefits of the implants. For augmentation, these outcomes included breast size change, satisfaction, and QoL measures. These outcomes were assessed before implantation and at 1, 2, and 3 years after surgery for those patients who still had their original implants and came back for follow-up visits.

**Primary Augmentation Patients:** For primary augmentation patients, 368 (67%) out of the original 552 patients were included in the analysis of cup size at 4 years. Of these 368 patients, 363 (99%) experienced at least one cup size increase. The average increase in circumferential chest size was 7.1 centimeters.

Mentor's satisfaction assessment was based on a single question of "Would the patient have this breast surgery again?" At 4 years, 394 (71%) of the 552 primary augmentation patients enrolled answered that question. Of these 394 patients, 385 (98%) stated to their surgeon that they would have the breast implant surgery again.

With regard to QoL measures at 4 years, an increase in self-esteem was noted for patients after primary breast augmentation on the Rosenberg Self-Esteem Scale. There was an improvement in the overall score of the Body Esteem Scale and the Sexual Attractiveness Subscale and the Chest Score of the Body Esteem Scale also increased. The SF-36 is a collection of scales assessing mental and physical health, and there was no improvement in the SF-36 after primary augmentation. The Tennessee Self-Concept Scale (TSCS) is a survey completed by the patient that evaluates how the patient sees herself and what she does, likes, and feels. There was no change in the overall score for the TSCS.

**Revision-Augmentation Patients:** For revision-augmentation patients, 97 (67%) out of the original 145 patients were included in the analysis of circumferential chest size at 4 years. For these 97 patients, the average increase in circumferential chest size was 2.6 centimeters.

Mentor's patient satisfaction was based on a single question of "Would the patient have this breast surgery again?" At 4 years, 97 (67%) of the 145 revision-augmentation patients enrolled answered that question. Of these 97 patients, 93 (96%) stated to their surgeon that they would have the breast implant surgery again.

With regard to QoL measures at 4 years, no change in self-esteem was noted following revision-augmentation surgery on the Rosenberg Self-Esteem Scale or the Body Esteem Scale, but the Chest Score of the Body Esteem Scale was significantly improved. The SF-36 is a collection of scales assessing mental and physical health, and there was no improvement in SF-36. The Tennessee Self-Concept Scale (TSCS) is a survey completed by the patient that evaluates how the patient sees herself and what she does, likes, and feels. There was no change in the overall TSCS score.

### **3.4. What Were the 4-Year Complication Rates in the Mentor Round Gel Core Study Augmentation Patients?**

The 4-year complication rates are shown from the most common to the least common in Table 1 (primary augmentation) and Table 2 (revision-augmentation) below. The rates reflect the percentage of augmentation patients who experienced the listed complication at least once within the first 4 years after implantation.

Some complications occurred more than once for some patients. The two most common complications experienced by primary augmentation patients within the first 4 years of implantation were reoperation (16%) and nipple sensation changes (12%).

**Table 1. 4-Year Complication Rates for the Mentor Round Gel Core Study Primary Augmentation Patients N=552 Patients**

Key Complications	%	CI
Reoperation	16.0	13.2, 19.4
Capsular Contracture Baker Grade III/IV	8.8	6.7, 11.6
Implant Removal with Replacement with Study Device	2.8	1.7, 4.6
Implant Removal with Replacement with Unknown Device	0	–
Implant Removal without Replacement	2.4	1.3, 4.1
Infection	1.5	0.7, 2.9
Rupture <sup>1</sup>	1.8	0.6, 5.5
Complications > 1% <sup>2</sup>	%	CI
Nipple Complications <sup>3</sup>	12.0	9.5, 15.1
Hypertrophic Scarring	6.8	5.0, 9.3
Other (Non-cosmetic) <sup>4</sup>	5.8	4.0, 8.3
Breast Mass	4.0	2.6, 6.1
Ptosis	3.0	1.8, 4.9
Breast Sensation Changes <sup>3</sup>	2.6	1.6, 4.4
Hematoma	2.5	1.5, 4.3
Breast Pain <sup>3</sup>	1.9	1.0, 3.5
Miscarriage	1.7	0.9, 3.3
Capsular Contracture Baker Grade II with Surgical Intervention	1.1	0.5, 2.4
External Injury Not Related to Breast Implants	1.0	0.4, 2.3

<sup>1</sup> There were 2 patients with suspected ruptures by MRI that were later determined to be no ruptures, and 1 patient who had suffered a hard fall with a confirmed rupture at explant in the augmentation group at 4 years.

<sup>2</sup> The following complications occurred at a rate less than 1%: asymmetry, breast lesion, car accident with seatbelts tender and burning, creasing of implants with palpable fold points, deep vein thrombosis, distortion of breast not related to capsular contraction, ecchymosis, excessive implant movements, excessive implant movement with muscle flexion, explanted due to right side being removed, fibrocystic changes, granuloma, implant malposition/displacement, implant palpability, implant riding high, inadequate milk supply, inflammation, inframammary fold too high, lactation difficulties, Mondor's disease, necrosis, new diagnosis of rheumatic disease, placement damage, positive nuclear antibodies negative for lupus, rash, right Montgomery gland infection, rippling, seroma, stillborn birth, suture reaction, thin spot in muscle wall, and wrinkling.

<sup>3</sup> Mild occurrences were excluded.

<sup>4</sup> Any complication other than ptosis, hypertrophic scarring, asymmetry, or wrinkling.

The two most common complications experienced by patients within the first 4 years of revision-augmentation surgery were reoperation (28%) and capsular contracture III/IV (20%). Notice that the rates for these two complications are higher than for primary augmentation. (For primary augmentation, reoperation was 16% and capsular contracture was 9%.)

**Table 2. 4-Year Complication Rates for the Mentor Round Gel Core Study Revision-Augmentation Patients N=145 Patients**

<b>Key Complications</b>	<b>%</b>
Reoperation	28.4
Capsular Contracture Baker Grade III/IV	19.9
Rupture <sup>1</sup>	7.3
Implant Removal with Replacement with Study Device	8.2
Implant Removal with Replacement with Unknown Device	0.9
Implant Removal without Replacement	5.0
Infection	1.4
<b>Complications &gt; 1%<sup>2</sup></b>	<b>%</b>
Nipple Complications <sup>3</sup>	12.3
Hypertrophic Scarring	8.6
Breast Mass	7.2
Other (Non-cosmetic) <sup>4</sup>	6.9
Capsular Contracture Baker Grade II with Surgical Intervention	5.1
Miscarriage	2.9
Hematoma	2.8
Ptosis	2.3
Breast Sensation Changes <sup>3</sup>	2.1
Seroma	2.1
Delayed Wound Healing	2.1
Breast Pain <sup>3</sup>	2.1
Wrinkling	1.4
Inflammation	1.4
Implant Malposition	1.4
Extrusion (intact)	1.4

<sup>1</sup> Three patients had a suspected rupture determined to be no rupture and 1 patient had a bilateral rupture confirmed at explant.

<sup>2</sup> The following complications occurred at a rate less than 1%: back and neck pain related to large implants, External injury not related to breast implant, lactation difficulties, muscle spasm, new diagnosis of rheumatic disease, palpability, shape change, siliconoma, skin lesion, surgical removal of ectopic pregnancy.

<sup>3</sup> Mild occurrences were excluded.

<sup>4</sup> Any complication other than ptosis, hypertrophic scarring, asymmetry, or wrinkling.

### 3.5. What Were the Main Reasons for Reoperation in the Mentor Round Gel Core Study Augmentation Patients?

There may be one or more reasons identified for having a reoperation (additional surgery after the primary or revision breast augmentation). Furthermore, there may be multiple surgical procedures (for example, implant removal with or without replacement, capsule procedures, incision and drainage, implant reposition, scar revision, etc.) performed during a reoperation. In Mentor's Round Gel Core Study, there were 181 additional surgical procedures performed in 113 reoperations involving 86 primary augmentation patients.

Table 3 below provides the main reason for each reoperation in primary augmentation patients following initial implantation that were performed through 4 years. The most common reason for reoperation through 4 years in primary augmentation patients was because of capsular contracture (36% of 113 reoperations).

**Table 3. Main Reasons for Reoperation in the Mentor Round Gel Core Study Primary Augmentation Patients through 4 Years**

Reason for Reoperation	n	% (of 113 Reoperations)
Capsular Contracture II, III, IV	41	36.3
Patient Request for Style/Size Change	15	13.3
Hematoma/Seroma	12	10.6
Scarring/Hypertrophic Scarring	12	10.6
Biopsy	9	8.0
Asymmetry	6	5.3
Ptosis (sagging)	4	3.5
Infection	3	2.7
Wrinkling	3	2.7
Delayed Wound Healing	2	1.8
Implant Malposition	2	1.8
Breast Pain	1	0.9
Extrusion of Intact Implant	1	0.9
Necrosis	1	0.9
Suspected Rupture <sup>1</sup>	1	0.9
<b>Total</b>	<b>113</b>	<b>100</b>

<sup>1</sup> The device was removed and found to be intact (not ruptured).

In Mentor's Round Gel Core Study, there were 110 additional surgical procedures performed in 60 reoperations involving 40 revision-augmentation patients. Table 4 below provides the main reason for each reoperation in revision-augmentation patients following initial implantation that were performed through 4 years. The most common reason for reoperation in revision-augmentation patients through 3 years was capsular contracture (40% of 60 reoperations).

**Table 4. Main Reasons for Reoperation in the Mentor Round Gel Core Study Revision-Augmentation Patients through 4 Years**

Reason for Reoperation	n	% (of 60 Reoperations)
Capsular Contracture II, III, IV	24	40.0
Patient Request for Style/Size Change	6	10.0
Biopsy	6	10.0
Hematoma/Seroma	5	8.4
Delayed Wound Healing	4	6.7
Scarring/Hypertrophic Scarring	3	5.0
Extrusion of Intact Implant	2	3.3
Implant Malposition	2	3.3
Asymmetry	2	3.3
Ptosis (sagging)	1	1.7
Infection	1	1.7
Wrinkling	1	1.7
Irritation/Inflammation	1	1.7
Rupture	1	1.7
Suspected Rupture <sup>1</sup>	1	1.7
<b>Total</b>	<b>60</b>	<b>100</b>

<sup>1</sup> The device was removed and found to be intact (not ruptured).

### 3.6. What Were the Main Reasons for Implant Removal in the Mentor Round Gel Core Study Augmentation Patients?

The main reasons for implant removal among primary augmentation patients in the Round Gel Core Study over the 4 years are shown in Table 5 below. There were 45 implants removed in 26 patients. Of these 45 implants, 24 (53%) were replaced. The most common reason for implant removal was patient request for style/size change (60% of the 45 implants removed).

**Table 5. Main Reasons for Implant Removal in the Mentor Round Gel Core Study Primary Augmentation Patients through 4 Years**

Reason for Removal	n	% (of 45 Explants)
Patient Request for Style/Size Change	27	60.0
Capsular Contracture III/IV	7	15.6
Breast Pain	2	4.4
Infection	2	4.4
Necrosis	2	4.4
Patient Request for Size Change	2	4.4

**Table 5. Continued on next page**

**Table 5. (Continued)**

Suspected Rupture <sup>1</sup>	1	2.2
Contralateral Explantation	1	2.2
Wrinkling	1	2.2
<b>Total</b>	<b>45</b>	<b>100</b>

<sup>1</sup> The device was removed and found to be intact (not ruptured).

The main reasons for implant removal among revision-augmentation patients in the Round Gel Core Study over the 4 years are shown in Table 6 below. There were 33 implants removed in 19 patients. Of these 33 implants, 18 (55%) were replaced. The most common reason for implant removal was patient request (36% of the 33 implants removed).

**Table 6. Main Reasons for Implant Removal in the Mentor Round Gel Core Study Revision-Augmentation Patients through 4 Years**

<b>Reason for Removal</b>	<b>n</b>	<b>% (of 33 Explants)</b>
Patient Request for Style/Size Change	12	36.4
Capsular Contracture III/IV	11	33.3
Rupture	2	6.1
Breast Mass	1	3.0
Explant to Match Other Side	1	3.0
Asymmetry	1	3.0
Capsular Contracture II	1	3.0
Extrusion of Intact Implant	1	3.0
Scarring/Hypertrophic Scarring	1	3.0
Infection	1	3.0
Suspected Rupture <sup>1</sup>	1	3.0
<b>Total</b>	<b>33</b>	<b>100</b>

<sup>1</sup> The device was removed and found to be intact (not ruptured).

### **3.7. What Were Other Clinical Data Findings in the Mentor Round Gel Core Study Augmentation Patients?**

Below is a summary of clinical findings from Mentor's Round Gel Core Study with regard to connective tissue disease (CTD), CTD signs and symptoms, cancer, lactation complications, reproduction complications, and suicide. These issues, along with others, are being further evaluated as part of a Mentor large post-approval study involving patients followed through 10 years.

### CTD Diagnoses

Three primary augmentation patients and one revision-augmentation patient in the Mentor Round Gel Core Study were reported to have a new diagnosis of CTD according to a rheumatologist. These diagnoses were Hashimoto's Thyroiditis at 2 years, rheumatoid arthritis at 3 years, hypothyroidism at 2 years, and a combined case of fibromyalgia and chronic fatigue syndrome at 3.5 years. It cannot be concluded that these CTD diagnoses were caused by the implants because there was no comparison group of similar women without implants.

### CTD Signs and Symptoms

Data on over 100 self-reported signs and symptoms, including about 50 self-reported rheumatological symptoms, were collected. Compared to before having the implants, significant increases were found for weakness, neck pain/stiffness, generalized aching, joint pain, frequent muscle pain, numbness of hands, combined fatigue, and combined pain in primary augmentation patients. No significant increases were found for the revision-augmentation patients. These increases were not found to be related to simply getting older over time. The Mentor Round Gel Core Study was not designed to evaluate cause and effect associations because there is no comparison group of women without implants, and because other contributing factors, such as medications and lifestyle/exercise, were not studied. Therefore, it cannot be determined whether these increases were due to the implants or not. However, you should be aware that you may experience an increase in these symptoms after receiving breast implants.

### Cancer

There were no primary augmentation patients with new diagnoses of breast cancer through 3 years in Mentor's Round Gel Core Study. As previous breast cancer was an exclusion criterion for primary augmentation patients, there were no reports of breast cancer recurrence in this indication. There were no reports of new diagnoses or recurrence of breast cancer in revision-augmentation patients. There were no reports of other cancers, such as brain, respiratory, or cervical/vulvar.

### Lactation Complications

Of the 15 primary augmentation patients who attempted to breastfeed following breast implantation in Mentor's Round Gel Core Study through 4 years, two (13%) experienced difficulty with breastfeeding and one (7%) reported inadequate milk. Of the 3 revision-augmentation patients who attempted to breastfeed after receiving breast implants, one (33%) had difficulty breastfeeding.

### Reproduction Complications

Nine (1.6%) of the primary augmentation patients, and 3 (2.1%) of the revision-augmentation patients in Mentor's Round Gel Core Study reported a miscarriage through 4 years.

### Suicide

There were no reports of suicide in either the primary augmentation or revision-augmentation indications in Mentor's Round Gel Core Study through 3 years.

## **4. Mentor CPG Core Study Results for Augmentation and Revision-Augmentation**

This section of this brochure summarizes the results of the Mentor CPG Core Study conducted on Mentor's CPG breast implants for primary augmentation and revision-augmentation. The Mentor CPG Core Study is the primary clinical study for this product. The results of the CPG Core Study give you useful information on the experience of other women with Mentor CPG breast implants. While the results cannot be used to predict your individual outcome, they can be used as a rough guide of what you may expect. Your own complications and benefits depend on many individual factors.

As a note, supplemental safety information was also obtained from the U.K. Sharpe/Collis Study and the literature to help assess long-term rupture rate and the consequences of rupture for Mentor's silicone gel-filled breast implants. The literature, which had the most available information on the consequences of rupture, was also used to assess other potential complications associated with silicone gel-filled breast implants. The key literature is referenced in throughout the Potential Breast Implant Complications section in this document.

### **4.1. Overview of the Mentor CPG Core Study**

The Mentor CPG Core Study is a 10-year study to assess safety and effectiveness in augmentation, reconstruction, and revision (revision-augmentation and revision-reconstruction) patients. Patient follow-up is at 10 weeks and annually through 10 years. Safety is being assessed via complications, such as implant rupture, capsular contracture, and reoperation. Benefit (effectiveness) is assessed by patient satisfaction and measures of quality of life (QoL).

The CPG Core Study consists of 955 patients, including 572 primary augmentation patients, 124 revision-augmentation patients, 191 primary reconstruction patients, and 68 revision-reconstruction patients. Of these, 252 primary augmentation patients, 56 revision-augmentation patients, 74 primary reconstruction patients, and 37 revision-reconstruction patients are in the MRI cohort, which means that they are assessed for silent rupture by MRI at years 1, 2, 4, 6, 8, and 10. The study is currently ongoing, with the results through 2 years reported in this brochure. Mentor will periodically update this brochure as more information becomes available. You should also ask your surgeon if he/she has any available updated clinical information.

Mentor's CPG Core Study results indicate that the risk of any complication or reoperation at some point through 3 years after implant surgery is 31% for primary augmentation patients and 39% for revision-augmentation patients. The information below provides more details about the complications and benefits you may experience.

Described below are the benefits and complications reported in the Mentor CPG Core Study for augmentation patients. The findings are described separately for primary augmentation and revision-augmentation patients.

### **4.2. What Were the 3-Year Follow-up Rates in the Mentor CPG Core Study Augmentation Patients?**

At the 3-year follow-up visit, data are reported for 84% of the eligible primary augmentation patients and 84% of the eligible revision-augmentation patients.

### 4.3. What Were the Benefits for the Mentor CPG Core Study Augmentation Patients?

The CPG Core Study measured a variety of outcomes that assessed the benefits of the implants. For augmentation, these outcomes included breast size change, satisfaction, and QoL measures. These outcomes were assessed before implantation and at 1, 2, and 3 years after surgery for those patients who still had their original implants and came back for follow-up visits.

**Primary Augmentation Patients:** For primary augmentation patients, 448 (78%) out of the original 572 patients were included in the analysis of cup size at 3 years. Of these 448 patients, 440 (98%) experienced at least one cup size increase. The average increase in circumferential chest size was 2.1 centimeters.

Mentor's satisfaction assessment was based on a single question of "Would the subject make the same decision to have this breast surgery?" At 3 years, 461 (81%) of the 572 patients enrolled answered that question. Of these 461 patients, 435 (94%) stated to their surgeon that they would have the breast implant surgery again.

In the primary augmentation cohort the General Health Perceptions score of the SF-36 worsened. There were no significant changes in any of the SF-36 scales. There was no change in the Health scale of the Body Esteem Scale, nor was there a change in the Weight Concern or Physical Condition scales, or the Overall scale. The Sexual and Chest scales were significantly improved. There was a significant increase from preoperative to 2 years postoperative in the Rosenberg Self-Esteem Scale. The Breast Evaluation Questionnaire (BEQ) is a questionnaire developed to assess breast satisfaction and quality of life outcomes among breast surgery patients. When asked how satisfied she was with the general appearance of her breasts, the large majority of primary augmentation patients (93%) said they "very satisfied or somewhat satisfied." For the three factors of 1) Comfort when not fully dressed, 2) Comfort when fully dressed, and 3) Satisfaction with breast attributes, there was a significant improvement in the Augmentation cohort for all three factors from preoperative to 2 years postoperative.

**Revision-Augmentation Patients:** For revision-augmentation patients, 94 (76%) out of the original patients were included in the analysis at 3 years. For these 94 patients, the average increase in circumferential chest size was 0.6 centimeters.

Mentor's patient satisfaction was based on a single question of "Would the subject make the same decision to have this breast surgery?" At 3 years, 98 (79%) of the 124 revision-augmentation patients enrolled answered that question. Of these 98 patients, 92 (94%) stated to their surgeon that they would have the breast implant surgery again.

With regard to body image/self-esteem measures at 2 years, the SF-36 Physical Functioning scale showed no significant results for treatment effect. In the Body Esteem Scale, the Chest scale increased. The other five scales showed no difference. For the Rosenberg scale there was no change. When asked in the BEQ how satisfied she was with the general appearance of her breasts, the large majority of revision augmentation patients (76%) said they "very satisfied or somewhat satisfied." The BEQ showed significant improvement in all three factors in these patients.

#### 4.4. What Were the 3-Year Complication Rates in the Mentor CPG Core Study Augmentation Patients?

The 3-year complication rates are shown in Table 7 (primary augmentation) and Table 8 revision-augmentation below. The rates reflect the percentage of augmentation patients who experienced the listed complication at least once within the first 3 years after implantation. Some complications occurred more than once for some patients. The two most common complications experienced by primary augmentation patients within the first 3 years of implantation were reoperation (13%) and ptosis (7%).

**Table 7. 3-Year Complication Rates for the Mentor CPG Core Study Primary Augmentation Patients  
N=572**

<b>Key Complications</b>	<b>%</b>
Reoperation	13.4
Capsular Contracture III / IV	1.3
Implant Removal with Replacement with Study Device	1.9
Implant Removal without Replacement	3.0
Infection	0.9
Implant Rotation	1.1
Rupture (MRI Cohort)	0
<b>Other complications occurring in <math>\geq</math> 1% of patients<sup>1</sup></b>	<b>%</b>
Ptosis	6.7
Nipple Sensation Changes <sup>2</sup>	3.8
Patient Requested Size Change	3.2
Breast Sensation Changes <sup>2</sup>	2.5
Hypertrophic Scarring	2.4
Scarring	2.2
Breast Pain <sup>2</sup>	2.2
Mass/Cyst	1.9
Wrinkling <sup>2</sup>	1.8
Position Dissatisfaction <sup>2</sup>	1.5
Capsular Contracture III	1.3
Breast Mass/Cyst Not Associated with Breast Implant	1.1
Hematoma	1.0

<sup>1</sup> The following complications occurred at a rate less than 1%: Baker II capsular contracture w/surgical intervention, Baker IV capsular contracture, bruising of breast, calcification, delayed wound healing, fluid oozing from wound, granuloma, implant movement, implant outline visible through skin, inframammary fullness, irritation/inflammation, lactation difficulties, loss of definition of inframammary fold, loss of upper fill, metastatic disease, miscarriage, new diagnosis of breast cancer, new diagnosis of rheumatic disease, nipple complication, palpability-implant, patient dissatisfied with aesthetic appearance of breast, patient dissatisfied with feel of implant, pectoralis movement with muscle contraction, rash, seroma, skin lesion, skin paresthesia, stitch complication, stretching of soft tissue, swelling (excessive), thickened capsule, tightness of skin over implant, and wound dehiscence.

<sup>2</sup> Mild occurrences were excluded.

The two most common complications experienced by patients within the first 3 years of revision-augmentation surgery were reoperation (15%) and implant removal without replacement (8%). Notice that the rates for reoperation are higher than for primary augmentation. (For primary augmentation, reoperation was 13% and implant removal without replacement was 3.0%.)

**Table 8. 3-Year Complication Rates for the Mentor CPG Core Study Revision-Augmentation Patients N=124 Patients**

<b>Key Complications</b>	<b>%</b>
Reoperation	15.3
Capsular Contracture III / IV	4.4
Implant Removal with Replacement with Study Device	4.4
Implant Removal without Replacement	7.6
Infection	0.8
Implant Rotation	1.8
Rupture (MRI Cohort)	1.1
<b>Other complications occurring in <math>\geq</math> 1% of patients<sup>1</sup></b>	<b>%</b>
Patient Requested Size Change	5.9
Nipple Sensation Changes <sup>2</sup>	5.8
Wrinkling <sup>2</sup>	5.0
Ptosis	4.7
Hypertrophic Scarring	3.5
Breast Mass/Cyst Not associated with Breast Implant	2.9
Capsular Contracture III	2.6
Implant Palpability <sup>2</sup>	2.6
Patient Dissatisfied with Feel of Implant	2.6
Wound Dehiscence	2.5
Breast Sensation Changes	2.2
Capsular Contracture IV	1.8
Position Dissatisfaction	1.8
Asymmetry <sup>2</sup>	1.7
Size Change (Physician Assessment Only)	1.7
Patient Dissatisfied with Aesthetic Appearance of Breast	1.7
Lack of Projection	1.2

<sup>1</sup> The following complications occurred at a rate less than 1%: Baker II capsular contracture w/surgical intervention, breast pain, implant movement, implant outline visible through skin, irritation/inflammation, new diagnosis of breast cancer, new diagnosis of rheumatic disease, seroma, and stitch complication.

<sup>2</sup> Mild occurrences were excluded.

#### 4.5. What Were the Main Reasons for the Mentor CPG Core Study Reoperation in Augmentation Patients?

There may be one or more reasons identified for having a reoperation (additional surgery after the primary or revision breast augmentation). Furthermore, there may be multiple surgical procedures (for example, implant removal with or without replacement, capsule procedures, incision and drainage, implant reposition, scar revision, etc.) performed during a reoperation. In Mentor's CPG Core Study, there were 183 additional surgical procedures performed in 89 reoperations involving 73 primary augmentation patients.

Table 9 below provides the main reason for each reoperation in primary augmentation patients following initial implantation that were performed through 3 years. The most common reason for reoperation through 3 years in primary augmentation patients was patient requested size change (14% of 89 reoperations).

**Table 9. Main Reasons for Reoperation in the Mentor CPG Core Study Primary Augmentation Patients through 3 Years**

Reason for Reoperation	n	% (of 89 Reoperations)
Patient Requested Size Change	12	13.5
Breast Mass/Cyst Not associated with Breast Implant	8	9.0
Ptosis	7	7.9
Hematoma	5	5.6
Hypertrophic Scarring/Scarring	5	5.6
Mass/Cyst	4	4.5
Position Dissatisfaction	4	4.5
Asymmetry	4	4.5
Calcification	4	4.5
Infection	3	3.4
Reason Unknown	3	3.4
Capsular Contracture III	3	3.4
Nipple Complications	2	2.2
Breast Pain	2	2.2
Patient Dissatisfied with Aesthetic Appearance of Breast	2	2.2
Implant Rotation	2	2.2
New Diagnosis of Breast Cancer	2	2.2
Scarring	2	2.2
Skin Lesion	2	2.2
Seroma	2	2.2
Wrinkling	2	2.2
Nipple – Unacceptably Low Sensitivity	1	1.1
Implant Movement	1	1.1
Delayed Wound Healing	1	1.1

**Table 9. Continued on next page**

**Table 9. (Continued)**

Granuloma	1	1.1
Irritation/Inflammation	1	1.1
Loss of Definition of Inframammary Fold	1	1.1
Capsular Contracture II with Surgical Intervention	1	1.1
Thickened Capsule	1	1.1
Wound Dehiscence	1	1.1
<b>Total</b>	<b>89</b>	<b>100</b>

In Mentor's CPG Core Study, there were 58 additional surgical procedures performed in 22 reoperations involving 18 revision-augmentation patients. Table 10 below provides the main reason for each reoperation in revision-augmentation patients following initial implantation that were performed through 3 years. The most common reason for reoperation in revision-augmentation patients through 3 years was wound dehiscence (18% of 22 reoperations).

**Table 10. Main Reasons for Reoperation in the Mentor CPG Core Study Revision-Augmentation Patients through 3 Years**

<b>Reason for Reoperation</b>	<b>n</b>	<b>% (of 22 Reoperations)</b>
Wound Dehiscence	4	18.2
Patient Requested Size Change	3	13.6
Position Dissatisfaction	2	9.1
Ptosis	2	9.1
Implant Movement	1	4.5
Asymmetry	1	4.5
Breast Pain	1	4.5
Delayed Wound Healing	1	4.5
Hypertrophic Scarring	1	4.5
New Diagnosis of Breast Cancer	1	4.5
Patient Dissatisfied with Appearance of Breast	1	4.5
Rupture	1	4.5
Size Change - Physician Assessment Only	1	4.5
Capsular Contracture III	1	4.5
Wrinkling	1	4.5
<b>Total</b>	<b>22</b>	<b>100</b>

#### 4.6. What Were the Main Reasons for Implant Removal in the Mentor CPG Core Study Augmentation Patients?

The main reasons for implant removal among primary augmentation patients in the CPG Core Study over the 3 years are shown in Table 11 below. There were 51 implants removed in 26 patients. Of these 51 implants, 19 (37%) were replaced. The most common reason for implant removal was patient requested size change (43% of the 51 implants removed).

**Table 11. Main Reasons for Implant Removal in the Mentor CPG Core Study Primary Reconstruction Patients through 3 Years**

Reason for Removal	n	% (of 51 Explants)
Patient Request for Size Change	22	43.1
Asymmetry	5	9.8
Position Dissatisfaction	4	7.8
Wrinkling	4	7.8
Patient Dissatisfied with Aesthetic Appearance of Breast	4	7.8
Capsular Contracture III/IV	3	5.9
Capsular Contracture II with Surgical Intervention	3	5.9
Implant Rotation	2	3.9
Ptosis	2	3.9
Reason Unknown	2	3.9
Total	51	100

The main reasons for implant removal among revision-augmentation patients in the CPG Core Study over the 3 years are shown in Table 12 below. There were 24 implants removed in 13 patients. Of these 24 implants, 9 (38%) were replaced. The most common reason for implant removal was patient request for size change (29% of the 24 implants removed).

**Table 12. Main Reasons for Implant Removal in the Mentor CPG Core Study Revision-Augmentation Patients through 3 Years**

Reason for Removal	n	% (of 24 Explants)
Patient Request for Size Change	7	29.2
Asymmetry	2	8.3
Patient Dissatisfied with Aesthetic Appearance of Breast	2	8.3
Wrinkling	2	8.3
Capsular Contracture III/IV	2	8.3
Ptosis	2	8.3
Reason Unknown	1	4.2
Breast Pain	1	4.2

**Table 12. Continued on next page**

**Table 12.** *(Continued)*

New Diagnosis of Breast Cancer	1	4.2
Wound Dehiscence	1	4.2
Position Dissatisfaction	1	4.2
Rupture	1	4.2
Size Change - Physician Assessment Only	1	4.2
Total	24	100

#### **4.7. What Were Other Clinical Data Findings in the Mentor CPG Core Study Augmentation Patients?**

Below is a summary of clinical findings from Mentor's CPG Core Study with regard to connective tissue disease (CTD), rheumatological signs and symptoms, cancer, neurological disease, neurological signs and symptoms, lactation complications, reproduction complications, and suicide. These issues, along with others, are being further evaluated as part of a Mentor large post-approval study of Round Gel implants involving patients followed through 10 years.

##### CTD Diagnoses

No primary augmentation patients and one revision-augmentation patient in the Mentor CPG Core Study were reported to have a new diagnosis of CTD according to a rheumatologist. This diagnosis was rheumatoid arthritis at 11 months. It cannot be concluded that this CTD diagnosis was caused by the implants because there was no comparison group of similar women without implants.

##### CTD Signs and Symptoms

Data on over 100 self-reported signs and symptoms, including about 50 self-reported rheumatological symptoms, were collected. Compared to before having the implants, statistically significant reductions were found in primary augmentation patients for heart murmur and cardiovascular. For revision-augmentation patients, significant increases were found for combined pain, muscle, and joint. The Mentor CPG Core Study was not designed to evaluate cause and effect associations because there is no comparison group of women without implants, and because other contributing factors, such as medications and lifestyle/exercise, were not studied. Therefore, it cannot be determined whether this increase was due to the implants or not, based on the CPG Core Study. However, you should be aware that you may experience an increase in these symptoms after receiving breast implants.

##### Cancer

There were two primary augmentation patients and one revision-augmentation patient with new diagnoses of breast cancer through 3 years in the Mentor CPG Core Study. As previous breast cancer was an exclusion criterion for primary augmentation patients, there were no reports of breast implant reoccurrence in this indication. There were no reports of other cancers, such as brain, respiratory, or cervical/vulvar.

### Neurological Disease, Signs, and Symptoms

In the Mentor CPG Core Study, data on patient-reported neurological disease and signs and symptoms were collected at each postoperative visit. Examples of neurological symptoms included weakness, numbness of feet, ringing in ears, and fatigue. Investigators also performed a physical exam at each postoperative visit.

Through 3 years, there were no reports of patients with neurological diseases in any indication. Compared to before implantation, significant increases were found for the combined category of pain in revision-augmentation patients through 3 years. No significant increases for neurological signs and symptoms were found in were found in primary augmentation patients.

### Lactation Complications

Three (16%) of the 19 primary augmentation patients who attempted to breastfeed following breast implantation in Mentor's CPG Core Study through 3 years experienced difficulty with breastfeeding. The two revision-augmentation patients who attempted to breastfeed after receiving breast implants had no difficulty breastfeeding.

### Reproduction Complications

There was 1 report of miscarriage for a primary augmentation patient. There were no reports of miscarriage in revision-augmentation patients in Mentor's CPG Core Study through 3 years.

### Suicide

There were no reports of suicide in primary augmentation or revision-augmentation patients in Mentor's CPG Core Study through 3 years.

## **5. Surgery Considerations for Receiving Breast Implants**

This section provides a discussion of surgical considerations for breast augmentation.

### **5.1. Surgical Considerations for Breast Augmentation**

#### **5.1.1. What Are the Alternatives to Breast Augmentation with Silicone Breast Implants?**

For primary augmentation patients, alternatives may include:

- Accept your breasts as they are and have no surgery.
- Wear a padded bra or external prostheses.
- Have mastopexy surgery (breast lift) without an implant.
- Have surgery with saline implants.

For revision-augmentation patients, alternatives may include:

- No revision
- Removal with or without replacement.

#### **5.1.2. Choosing a Surgeon**

When choosing a surgeon who is experienced with breast augmentation, you should know the answers to the following types of questions:

- Has he/she received professional training, including techniques specific to plastic surgery and breast implantation?
- How many breast augmentation implantation procedures does he/she perform per year?
- How many years has he/she performed breast augmentation procedures?
- Is he/she board certified, and if so, with which board?
- What is the most common complication he/she encounters with breast augmentation?
- What is his/her reoperation rate with breast augmentation, and what is the most common type of reoperation he/she performs?
- Can he/she perform this surgery in a hospital, as well as in the surgeon's independent surgery center? (Note that hospitals require the demonstration of evidence of appropriate training in specific procedures before allowing surgeons to operate in their facilities.)

### **5.1.3. Implant Shape and Size**

Depending on the desired shape you wish to achieve, you and your surgeon have implants with three different round profiles, or styles, from which to choose. Generally, the larger you want your cup size, the larger the breast implant the surgeon will consider (measured in cubic centimeters, or cc's), not in cup sizes, because this depends on the size and shape of the individual woman's chest.

Your surgeon will also evaluate your existing breast and skin tissue to determine if you have enough to cover the breast implant you are considering, or, in some cases such as after pregnancy, too much extra skin. If you desire a breast implant size that is too large for your tissue, the surgeon may warn you that breast implant edges may be visible or palpable postoperatively. Also, excessively large breast implants may speed up the effects of gravity on the breast and can result in droop or sag at an earlier age. One report indicates that larger sized implants (greater than 350 cc) may be too large for many women, increasing the risk of developing complications such as implant extrusion, hematoma, infection, palpable implant folds, and visible skin wrinkling requiring surgical intervention to correct these complications.<sup>79</sup>

### **5.1.4. Surface Texturing**

Some studies suggest that surface texturing reduces the chance of severe capsular contracture,<sup>80</sup> while other studies do not.<sup>81,82</sup> Mentor's Round Gel Core Study did not show a difference in the likelihood of developing capsular contracture with textured implants compared to smooth-surfaced implants. Mentor's CPG implants are only available with a textured surface.

A textured implant may require a larger incision because the rougher textured surface makes it harder to place into the pocket without undue stress, which might damage the implant or decrease its durability.

The risk of BIA-ALCL is higher for textured surface breast implants versus smooth surface breast implants.

### 5.1.5. Materials Used in the Manufacture of Breast Implants

Below is a summary of materials found in MemoryGel Breast Implants and CPG Breast Implants


#### MemoryGel and CPG Breast Implant Device Materials:

Device Materials	Implant Component
Dimethyl Silicone Elastomer Dispersion	Shell
Diphenyl Silicone Elastomer Dispersion	Shell
Silicone Elastomers	Shell
Silicone Gel	Gel fill
Platinum (crosslinking catalyst)	Shell and fill

The potential toxicity of the materials listed in the above table has been evaluated with both toxicity testing and risk assessments to assess the exposure levels in comparison to the amount determined to likely be safe. However, individual responses to substances may vary, and all reactions cannot be predicted. Further information on components that can be extracted from breast implants can be found on FDA's website ([https://www.accessdata.fda.gov/cdrh\\_docs/pdf3/p030053b.pdf](https://www.accessdata.fda.gov/cdrh_docs/pdf3/p030053b.pdf)).

### 5.1.6. Implant Placement

The breast implant can be placed either partially under the pectoralis major muscle (submuscular) or on top of the muscle and under the breast glands (subglandular). You should discuss with your surgeon the advantages and disadvantages of the implant placement selected for you. Each implant location has benefits and negative features of which you need to be aware, as described in Table 13 below.

**Table 13. Comparison between Submuscular versus Subglandular Placement**


Submuscular Placement	Subglandular Placement
Surgery may be longer	Surgery may be shorter
Recovery may be longer	Recovery may be shorter
May be more painful	May be less painful
Reoperation may be more difficult	May provide easier access for reoperation
Less visible and palpable implants	More visible and palpable implants
Less likelihood of capsular contracture <sup>83</sup>	Greater likelihood of capsular contracture <sup>84,85</sup>
Easier imaging during mammography exam	More difficulty imaging during mammography exam
May be preferable if you have thin or weakened breast tissue	May not be recommended if you have thin or weakened breast tissue.

### 5.1.7. Incision Sites

You should discuss with your surgeon the pros and cons for the incision site specifically recommended for you.

The incision size will be larger than for a saline breast augmentation. There are 3 common incision sites: under the arm (axillary), around the nipple (periareolar), or within the breast fold (inframammary).

- **Periareolar** – This incision is typically more concealed, but since it also involves cutting through the breast tissue, it is associated with a higher likelihood of breastfeeding difficulties, as compared to the other incision sites.<sup>86</sup> Cutting through the tissue may increase the chance that there will be a change in breast and/or nipple sensation.
- **Inframammary** – This incision is generally less concealed than periareolar and associated with less breastfeeding difficulties than the periareolar incision site. It is also the most commonly used incision site at the present time, and is felt to give the best access to and control of the breast implant pocket.
- **Axillary** – This incision is less concealed than periareolar and associated with less breastfeeding difficulties than the periareolar incision site. If the incision is made under the arm, the surgeon may use a probe fitted with a miniature camera, along with minimally invasive (very small) instruments, to create a “pocket” for the breast implant. This approach is more difficult, and may increase the risk of damage to, and unexpected location of, the implant.
- **Umbilical** (belly button) – This incision site has not been studied in Mentor’s Round Gel or CPG Core Studies and should not be used for a wide variety of reasons, including potential damage to the implant shell.

### 5.1.8. Additional Procedures at the Time of Breast Augmentation

Your surgeon will examine your breasts and help you make decisions to obtain the best result in your individual situation. In some cases, particularly after pregnancy or significant weight loss, implants alone may not address all of the issues, such as sagging or extra skin, affecting your breasts. This is particularly true when there is extra skin remaining from when the breasts were engorged with milk, or when you might have been carrying more weight.

In these situations, your surgeon may recommend a breast lift (mastopexy) to remove some of the extra skin, or to lift the breasts, at the time of implant placement. Mastopexy involves removing a strip of skin from under the breast or around the nipple to lift the nipple and breast location and tighten the skin over the breast. Your surgeon will discuss the potential risks, and the location of the additional scars which might be required to lift your breasts or to remove the extra skin.

### 5.1.9. Palpability

Implants may be more palpable or noticeable if there is an insufficient amount of skin/tissue available to cover the implant and/or when the implant is placed subglandularly.

### **5.1.10. Surgical Setting and Anesthesia**

Augmentation surgery is usually performed on an outpatient basis, in a specialized operating room that may be located in a hospital, a surgery center, or surgical suite in the surgeon's office. General anesthesia is commonly used, and local anesthesia with sedation is also an option. You should be sure to check with your surgeon and with the facility where the surgery will take place, to become aware of the tests, presurgical examinations, and length of time you need to be without food or your routine medications prior to the surgical procedure.

### **5.1.11. Postoperative Care**

You will probably feel somewhat tired and sore for several days following the operation, and your breasts may remain swollen and sensitive to physical contact for a month or longer. You may also experience a feeling of tightness in the breast area as your skin adjusts to your new breast size. The feeling in the breasts and nipple area also may be diminished during this time of swelling and immediate post-surgery recovery. Other possible complications have been described above.

Postoperative care depends on each patient's situation, may involve the use of a special postoperative bra, compression bandage, or jog bra for extra support and positioning while you heal. Some surgeons may not want you to wear a bra at all for a period of time following the surgery. At your surgeon's recommendation, you will most likely be able to return to work within a few days, although for at least a couple of weeks you should avoid any strenuous activities that could raise your pulse and blood pressure or require strenuous use of your arms and chest. Your surgeon may also recommend breast massage exercises.

Note: If you experience fever, do not feel well, or see noticeable swelling and/or redness or drainage in your implanted breast(s), you should contact your surgeon immediately.

## **5.2. Other Factors to Consider in Revision-Augmentation Surgery**

Some revision surgeries require removal of an intact implant (for example, capsulotomy and pocket adjustments), while others do not require removal of the implant. Any device that has been removed during revision surgery should not be reimplanted. Mentor breast implants are "for single use only."

## **6. Follow-Up Examinations**

### **6.1. Breast Self-Examinations**

You should perform a breast self-examination monthly. This may be more difficult with an implant in place. In order to do this effectively, you should ask your surgeon to help you tell the difference between the implant and your breast tissue. Care should be taken not to squeeze the implant excessively. Any new lumps may be evaluated with a biopsy, as appropriate. If a biopsy is performed, care must be taken to avoid injuring the implant.

## **6.2. Recommended Imaging Schedule for Implant RUPTURE Surveillance**

Imaging is used for breast cancer screening and, for women with silicone gel implants, for breast implant rupture surveillance. Breast cancer screening can be routine or individualized depending on patient and family history and should be directed by your physician. The guidelines for surveillance of breast implant rupture are as follows:

It is recommended that you have periodic imaging (ultrasound or magnetic resonance imaging (MRI)) of your silicone gel-filled breast implants to screen for implant rupture regardless of whether your implants are for cosmetic augmentation or reconstruction. These recommendations do not replace other additional imaging that may be required depending on your medical history or circumstances (i.e., screening mammography for breast cancer).

Even if you have no symptoms, you should have your first ultrasound or MRI at 5-6 years after your initial implant surgery and then every 2-3 years thereafter. If you have symptoms at any time or uncertain ultrasound results for breast implant rupture, an MRI is recommended.

## **6.3. Mammography**

The current recommendations for getting screening/preoperative mammograms are no different for women with breast implants than for those without implants. Mammography exams should be interpreted by radiologists experienced in the evaluation of women with breast implants. It is essential that you tell your mammography technologist that you have an implant before the procedure. You should request a diagnostic mammogram, rather than a screening mammogram, because more pictures are taken with diagnostic mammography. The technologist can use special techniques to reduce the possibility of rupture and to get the best possible views of the breast tissue. More information on mammography is provided in Section 1.4.

## **7. The Types of Silicone Gel-Filled Breast Implants Available from Mentor**

Mentor's silicone gel-filled breast implants, referred to as MemoryGel products, come in a variety of profiles and sizes. All currently available MemoryGel breast implants have either a textured shell or smooth surface shell. MemoryGel Xtra and MemoryGel BOOST Breast Implants are available in a smooth shell only.

Mentor's CPG breast implants also come in a variety of profiles and sizes. All currently available CPG breast implants have a textured shell surface. Be sure to familiarize yourself with the different features of breast implants and to discuss the best type(s) of implants for you with your surgeon.

## **8. Product Replacement Policy and Limited Warranties**

The following is a description of the assistance available from the Mentor Lifetime Product Replacement Policy and the Mentor Advantage Limited Warranty.

### **Mentor's Free Lifetime Product Replacement Policy**

- Automatically applies to all recipients of Mentor breast implant products.
- Provides that regardless of the age of the implant, when confirmed deflation or rupture occurs, you are eligible for 1 to 2 no charge replacement breast implant products of any size in a similar style.

The **Mentor Advantage** is free of charge to all patients who are implanted with Mentor saline-filled breast or silicone gel-filled implant products. When the limited warranty applies, Mentor provides the following:

- Lifetime product replacement policy\*
- 10 years and up to \$1200(CAD) financial assistance for operating room, anesthesia, and surgical charges not covered by insurance\*\*
- Free contralateral (opposite side) implant replacement upon surgeon request
- Non-cancellable terms.

With the Mentor Advantage Warranty, it is important for you to also maintain your own records to ensure validation of your enrollment, as it is possible your surgeon may not retain your records for the entire duration of the limited warranty.

### **Products Covered**

The Mentor Advantage coverage applies to all Mentor breast implants that are implanted in Canada after May 1, 2005<sup>†</sup>, provided implants have been:

- Implanted in accordance with the Mentor package insert, current to the date of implantation, and other notifications or instructions published by Mentor; and
- Used by appropriately qualified, licensed surgeons, in accordance with accepted surgical procedures.

### **Events Covered**

The Mentor Advantage coverage applies to the following:

- Deflation due to crease fold failure, patient trauma, or unknown cause
- Loss of valve integrity
- Other loss-of-shell-integrity events also may be covered by this program. Mentor reserves the right to determine if specific, additional, events should be covered.

### **Events Not Covered**

The Mentor Advantage does not apply to the following:

- Removal of intact implants due to capsular contracture, wrinkling, or rippling.
- Loss of implant shell integrity resulting from reoperative procedures, open capsulotomy, or closed compression capsulotomy procedures.
- Removal of intact implants for size alteration.

\* Lifetime Product Replacement Policy: Mentor will provide replacement Mentor product of any size in the same or similar style as the originally implanted product free of charge for the lifetime of the patient. Upon surgeon's request, a different implant style maybe selected (subject to a charge of the difference between product list prices). See the Mentor Advantage Limited Warranty for eligibility and program details.

\*\* Operating room and anesthesia charges to be given payment priority. In order to qualify for financial assistance, you will need to sign a Release form.

† For breast implants implanted prior to this date, contact Mentor Worldwide for information regarding any applicable warranty terms.

## Filing for Financial Assistance

- To file a Mentor Advantage claim for product replacement and/or financial assistance, the surgeon must contact Mentor Medical Systems Canada prior to replacement surgery.
- For financial assistance claims, a patient-specific Release form will be generated that you must sign and return.
- For either replacement or financial assistance claims, the surgeon must send the explanted, decontaminated Mentor breast implant(s) within six months of the date of explantation (implant removal) to:  
Mentor, a unit of Johnson & Johnson Medical Products, a division of Johnson & Johnson Inc.  
200 Whitehall Drive  
Markham, ON Canada  
L3R 0T5
- Upon receipt, review and approval of the completed claim, including receipt of the explanted product and your completion of a full general release, financial assistance will be issued.

This is a summary of the coverage of the Mentor Advantage Limited Warranty. It is an overview only and not a complete statement of the program. A copy of the complete Mentor Advantage Limited Warranty for saline-filled and silicone gel-filled breast implants may be obtained by writing or calling:

Mentor, a unit of Johnson & Johnson Medical Products, a division of Johnson & Johnson Inc.  
200 Whitehall Drive  
Markham, ON  
Canada L3R 0T5  
1-800-668-6069

A copy of the complete programs may also be obtained from your surgeon or by going to [www.mentorwwllc.com](http://www.mentorwwllc.com).

**THIS IS A LIMITED WARRANTY ONLY AND IS SUBJECT TO THE TERMS AND CONDITIONS SET FORTH AND EXPLAINED IN THE APPLICABLE MENTOR LIMITED WARRANTIES. ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, BY OPERATION OF LAW OR OTHERWISE, INCLUDING BUT NOT LIMITED TO, IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS ARE EXCLUDED.**

Mentor reserves the right to cancel, change, or modify the terms of the Mentor Advantage coverage. Any such cancellation, change, or modification will not affect the currently stated terms of the Mentor Advantage coverage for those already enrolled.

## **9. Other Sources of Additional Information**

Upon request, you will be provided with a copy of the package insert (Directions for Use). You can request a copy from your surgeon or from Mentor. The package insert has many undefined medical and technical terms since it contains information directed only to the surgeon.

If you should decide to get breast implants, you will be given a Patient Implant Identification card with the style and serial number of your breast implant(s). This will be given to you right after your surgery. It is important that you keep a copy of this card because you may need to see that information at a later date. A mammography center form and a primary physician form have been provided at the end of this brochure. These forms are available for you to record the information about the breast implants you have received and bring it with you to your next mammography center or primary physician appointments.

For additional information or questions about Mentor breast implants, please call 1-800-MENTOR8. Other resources for information about breast implants are provided below.

### **Mentor Worldwide LLC**

1-800-525-0245

[www.mentorwllc.com](http://www.mentorwllc.com) [www.loveyourlook.com](http://www.loveyourlook.com)

### **Institute of Medicine Report on the Safety of Silicone Implants**

[www.nap.edu/catalog/9618.html](http://www.nap.edu/catalog/9618.html)

### **Health Canada**

[http://www.hc-sc.gc.ca/iyh-vsv/med/implants\\_e.html](http://www.hc-sc.gc.ca/iyh-vsv/med/implants_e.html)

You can find important information in the Health Canada “It’s Your Health - Breast Implants” handbook, which is available through the website provided above.

To report any problems or adverse effects you have with your breast implants to Health Canada:

- Call toll-free at 1-800-267-9675
- Fill out the Consumer Medical Device Report Form to report a medical device problem

### **Canadian Society of Plastic Surgeons**

<http://www.plasticsurgery.ca>

## **10. PATIENT CHECKLIST – TO BE COMPLETED PRIOR TO SURGERY**

### **10.1. Introduction:**

Below is a checklist of the risks associated with breast implants and breast implant surgery. The checklist also includes other important information, such as health insurance coverage, for you to consider. Please take time and review each section of the checklist. It may be helpful to place your initials at the end of each section to indicate your understanding of the information presented or, if there are sections that you are unsure about, write down your questions and discuss them with your surgeon before deciding to have breast implant surgery.

This checklist can be used as a tool to inform you of the described risks of breast implants.

When you place your signature at the end of the checklist, you are confirming that you have reviewed each section, have had your questions addressed and understand all the information presented. You can ask your surgeon to sign the checklist and to make it part of your permanent medical record.

It is important to remember that the lifespan of breast implants varies by person and cannot be predicted. That means everyone with breast implants may need additional surgeries, but no one can predict when. The longer your implants are implanted, the greater the chances are that you will develop complications, some of which will require more surgery.

Mentor's Canadian patient website, [breastimplantsbymentor.net/en-CA](http://breastimplantsbymentor.net/en-CA), includes a Safety Information Section. The latest version of Mentor's Patient Educational Brochures, physician labeling (Product Insert Data Sheet) and updates on important safety information are available in this section. You should check this website periodically to stay up to date on any new safety information posted.

Choosing the right surgeon for your breast implant surgery is very important. You should ask your physician about his/her experience, medical degree, specialty of training and credentials so that you are assured your surgery will be performed by an appropriately trained surgeon.

## 10.2. Are you eligible for breast implant surgery?

There are some conditions that would prevent you from undergoing breast implant surgery or put you at higher risk for a poor surgical outcome.

You are not able to receive breast implants if you have an infection, cancer or pre-cancer of your breast tissue that has not been adequately treated, or if you are pregnant or nursing.

The risk of complications due to breast implant surgery may be higher if you are a smoker, have a weakened immune system, have undergone or plan to undergo chemotherapy or radiation therapy, have a condition that interferes with wound healing or blood clotting, or if you have a reduced blood supply to the breast tissue. The risk may also increase if you have an autoimmune disease or if you suffer from clinical depression or other mental disorders.

✓ **Place initials here if you understand the information presented in the section:** \_\_\_\_\_

➤ **Jot down your questions here if you do not understand this section and follow up with your surgeon**

## 10.3. Regular follow up and imaging of your breast implants

Even if you do not experience issues with your breast implants, you should have your breasts examined annually by your plastic surgeon or a physician that is qualified to examine breasts with implants.

Breast implants may make it more difficult to monitor for breast cancer during a routine mammogram and additional images may need to be taken. Mammogram may cause your implants to rupture or leak. If you have implants following breast reconstruction, mammograms are not recommended so these concerns do not apply.

Breast implants may rupture or leak over time. If you have gel-filled (silicone-filled) breast implants, you may not know that rupture has occurred. This is called “silent rupture.” because this can occur without symptoms or visible changes. Gel-filled breast implants should be checked for silent rupture beginning 5 years after your surgery and every 2 years after that. Your implants will need to be checked using ultrasound or Magnetic Resonance Imaging (MRI). These tests may not be covered by insurance. If rupture is detected, your implants should be removed.

Breast implants contain various materials including silicone polymers, platinum and other trace metals. Studies have determined that these materials are present at a safe level. You should make your surgeon aware of any sensitivities that you may have to specific substances.

If your gel-filled breast implant leaks or ruptures, it is possible that silicone from your implant may travel to nearby tissues (e.g., lymph nodes under the arm) and organs (e.g., liver, lungs) where it may not be possible to remove. The health consequences of released (migrated) medical grade silicone are not fully understood.

- ✓ **Place initials here if you understand the information presented in the section: \_\_\_\_\_**
- **Jot down your questions here if you do not understand this section and follow up with your surgeon**

#### **10.4. Insurance may not cover your breast implant surgery.**

Health insurance does not cover the cost of breast implants that are used to change the size or shape of your breasts for cosmetic purposes and it may not cover treatments for complications that may develop in the future.

- ✓ **Place initials here if you understand the information presented in the section: \_\_\_\_\_**
- **Jot down your questions here if you do not understand this section and follow up with your surgeon**

#### **10.5. Other options for augmenting your breasts**

If you are considering breast implants to enhance your appearance, you should be aware that there are other options for augmenting your breasts such as fat transplantation or the use of a padded bra.

- ✓ **Place initials here if you understand the information presented in the section: \_\_\_\_\_**
- **Jot down your questions here if you do not understand this section and follow up with your surgeon**

#### **10.6. Risks associated with all breast surgery**

There are risks that can occur with any surgery involving the breast. Some of the more commonly occurring complications include:

- Hematoma: bleeding at the surgical site after wound closure
- Seroma: excess fluid collection at the surgical site after wound closure
- Infection: bacterial contamination of the surgical site during healing

- Loss of breast skin or nipple sensation and/or erection due to nerve damage
- Changes in your breast shape may occur as your body heals after surgery and as your body ages
- Insufficient breast milk production: Some women with breast augmentation have reduced lactation capacity.
- Breast pain
- Asymmetry
- Swelling
- Scarring
- Tissue death of breast skin or nipple
- Anesthetic complications

✓ **Place initials here if you understand the information presented in the section:** \_\_\_\_\_

➤ **Jot down your questions here if you do not understand this section and follow up with your surgeon**

### 10.7. Risk associated with surgically placed breast implants

Complications are known to be associated with breast implants. Some can occur soon after surgery while others can occur years after implantation. Some of the more common risks known to occur are discussed below. You can find more detailed information, including the percentage of women in our clinical study who had each type of complication, in Section 3 and Section 4 of this brochure. You should ask your surgeon to review the clinical data if you have questions.

- **Additional Surgeries:** Breast implants are not lifetime devices, and you may need additional surgery to replace or remove your implants. You may also need additional surgery to address known complications that can develop over time (such as those discussed below) or to address possible dissatisfaction with your postoperative appearance, symmetry or size.
- **Rupture:** The shell of your implant may tear or leak over time. It is recommended that an implant known to be ruptured be removed. If you have a saline-filled breast implant(s), you will know if a rupture occurs because the implant loses volume and deflates as the body absorbs free saline. If you have a gel-filled breast implant(s), you may not experience volume loss with rupture. This is why routine screening for implant rupture with imaging is important.
- **Permanent Changes to Breast Tissue:** Breast implantation for augmentation can result in changes to the original breast skin envelope, glandular volume and nipple position that may be permanent despite future implant removal. In some cases, breast tissue atrophy can occur.

- **Capsular Contracture:** Scar tissue naturally forms around your implant. In some cases, the scar tissue may become too thick and organized and squeeze the implant. This can be painful and may give your breasts an unnatural appearance or firmness.
- **Implant Infection:** An infection may develop in the tissue around your implants and, in some cases, this can only be resolved by implant removal.
- **Implant Exposure or Implant Extrusion:** This occurs when the breast implant becomes exposed, usually due to the incision reopening at the insertion site.
- **Wrinkling/Rippling** can be visible from the skin, the capsule or the implant
- **Visibility of the Implant Edges Through the Breast Tissues**
- **Incorrect positioning of the Implant** either immediately or years after surgery

One complication, while it does not occur frequently, that you should be aware of is **Breast Implant Associated Anaplastic Large Cell Lymphoma (BIA-ALCL)**.

- BIA-ALCL is a cancer of the immune system. It occurs more commonly in patients with textured implants than smooth implants, although rates are not well defined. BIA-ALCL typically takes several years to develop after implantation but some cases have been reported as early as one year after implantation. Symptoms to be aware of include breast swelling, breast masses and lymph node enlargement. Some women with advanced disease at diagnosis have died from BIA-ALCL.
- You can also find more detailed information about BIA-ALCL can be found in Section 2 in of this brochure, on the FDA website at [www.fda.gov/medical-devices/breast-implants/questions-and-answers-about-breast-implant-associated-anaplastic-large-cell-lymphoma-bia-alcl](http://www.fda.gov/medical-devices/breast-implants/questions-and-answers-about-breast-implant-associated-anaplastic-large-cell-lymphoma-bia-alcl) and the Health Canada website at [www.canada.ca/en/health-canada/services/drugs-medical-devices/breast-implants](http://www.canada.ca/en/health-canada/services/drugs-medical-devices/breast-implants)

Some women with implants may have experienced health problems such as **Connective Tissue Diseases** (e.g., lupus, rheumatoid arthritis, scleroderma). There is currently insufficient evidence to support an association between implants and those diagnoses.

Some women receiving breast implants have reported a variety of **Systemic Symptoms (also referred to as Breast Implant Illness)** such as joint pain, muscle aches, confusion, chronic fatigue and other rheumatology symptoms. Some women have also reported relief of these symptoms, at least for some time, after their implants were removed. Individual patient risk for developing these symptoms has not been well established.

✓ **Place initials here if you understand the information presented in the section:** \_\_\_\_\_

➤ **Jot down your questions here if you do not understand this section and follow up with your surgeon**

**10.8. Acknowledgment of Understanding****Patient:**

I acknowledge that I have received and read the Patient Education Brochure for the implant that will be implanted as part of my breast augmentation surgery and understand the benefits and risks of breast implants for me, given my specific health conditions. I understand the alternatives to breast implant-based breast augmentation. Any questions related to this material have been adequately addressed with the health care professionals involved in my care.

---

Patient signature and date

**Physician:**

I acknowledge that I have discussed with my patient the benefits and risks of breast implants as described in this Patient Education Brochure and in this checklist. I have also explained the benefits and risks of the alternatives. I have addressed questions and feel the patient has a thorough understanding of the risks discussed in this brochure.

---

Physician's signature and date

## **11. MAMMOGRAPHY CENTER FORM**

### **For Patients:**

Remember to take this information about your breast implant surgery to your Mammography Center Appointments.

### **Information For The Mammography Center:**

Please update my patient file to reflect the presence of Mentor™ MemoryGel™, Mentor™ MemoryGel™ Xtra or Mentor™ MemoryGel™ CPG Silicone Gel-Filled Breast Implants. Since examination of women with breast implants requires more time, please allow additional time when scheduling my next mammogram and alert the physician and technologists performing the exam about the presence of my implants.

Health Canada has approved Mentor™ MemoryGel™, Mentor™ MemoryGel™ Xtra and Mentor™ MemoryGel™ CPG Silicone Gel-Filled Breast Implants for use in augmentation, reconstruction and revision surgery. As part of a woman's healthcare network, it is important that you are aware of the latest information on the safety of silicone gel-filled breast implants. For additional information or questions about Mentor breast implants, please call 1-800-MENTOR8. Other resources for information about breast implants are provided below.

### **Mentor Worldwide LLC**

1-800-525-0245

[www.mentorwwllc.com](http://www.mentorwwllc.com)

### **Health Canada**

[http://www.hc-sc.gc.ca/iyh-vsv/med/implants\\_e.html](http://www.hc-sc.gc.ca/iyh-vsv/med/implants_e.html)

### **Canadian Society of Plastic Surgeons**

<http://www.plasticsurgery.ca>

### **Institute of Medicine Report on the Safety of Silicone Implants**

[www.nap.edu/catalog/9618.html](http://www.nap.edu/catalog/9618.html)

**CANADIAN  
MAMMOGRAPHY CENTER FORM**

**FOR PATIENTS:**

If you should decide to get breast implants, you will be given a Patient Implant Identification Card with the style and serial number of your breast implant(s). This will be given to you right after your surgery. Bring this card with you to your Mammography Center or please record the catalog and serial number(s) exactly as they appear on your MENTOR Patient Implant Identification Card(s) before giving this page to your Mammography Center.

**RECORD RIGHT Implant Catalog Number and Serial Number below :**

**Catalog Number:** \_\_\_\_\_

**Serial Number:** \_\_\_\_\_

**RECORD LEFT Implant Catalog Number and Serial Number below :**

**Catalog Number:** \_\_\_\_\_

**Serial Number:** \_\_\_\_\_

## **12. PRIMARY PHYSICIAN FORM**

### **For Patients:**

Remember to provide this information about your implant surgery to your Primary Care Physician.

### **Information For Your Primary Care Physician:**

Your patient has been implanted with Mentor™ MemoryGel™, Mentor™ MemoryGel™ Xtra or Mentor™ MemoryGel CPG Silicone Gel-Filled Breast Implants. Due to the potential risks associated with silicone gel-filled breast implants, it is important that you include this information in her chart. To ensure your patient receives the care she needs, when appropriate please alert other physicians about the presence of her implants.

Health Canada has approved Mentor™ MemoryGel™, Mentor™ MemoryGel™ Xtra and Mentor™ MemoryGel CPG Silicone Gel-Filled Breast Implants for use in augmentation, reconstruction and revision surgery. As part of a woman's healthcare network, it is important that you are aware of the latest information on the safety of silicone gel-filled breast implants. For additional information or questions about Mentor breast implants, please call 1-800-MENTOR8. Other resources on breast implants are listed below.

### **Mentor Worldwide LLC**

1-800-525-0245

[www.mentorwwllc.com](http://www.mentorwwllc.com)

### **Health Canada**

[http://www.hc-sc.gc.ca/iyh-vsv/med/implants\\_e.html](http://www.hc-sc.gc.ca/iyh-vsv/med/implants_e.html)

### **Canadian Society of Plastic Surgeons**

<http://www.plasticsurgery.ca>

### **Institute of Medicine Report on the Safety of Silicone Implants**

[www.nap.edu/catalog/9618.html](http://www.nap.edu/catalog/9618.html)

**CANADIAN  
PRIMARY CARE PHYSICIAN FORM**

**FOR PATIENTS:**

If you should decide to get breast implants, you will be given a Patient Implant Identification Card with the style and serial number of your breast implant(s). This will be given to you right after your surgery. Bring this card with you to your Primary Care Physician or please record the catalog and serial number(s) exactly as they appear on your MENTOR Patient Implant Identification Card(s) before giving this page to your Primary Care Physician.

**RECORD RIGHT Implant Catalog Number and Serial Number below :**

**Catalog Number:** \_\_\_\_\_

**Serial Number:** \_\_\_\_\_

**RECORD LEFT Implant Catalog Number and Serial Number below :**

**Catalog Number:** \_\_\_\_\_

**Serial Number:** \_\_\_\_\_

## References

1. Bondurant, S., V.L. Ernster and R. Herdman, Eds. 2000. Safety of silicone breast implants. Committee on the Safety of Silicone Breast Implants, Division of Health Promotion and Disease Prevention, Institute of Medicine. Washington, D.C.: National Academy Press.
2. Hölmich, L.R., et al. 2004. Untreated silicone breast implant rupture. *Plast. Reconstr. Surg.* 114:204-214.
3. Hölmich, L.R., et al. 2003a. Incidence of silicone breast implant rupture. *Arch Surg.* 138:801-806.
4. Holmich, L.R., et al. 2001. Prevalence of silicone breast implant rupture among Danish women. *Plast. Reconstr. Surg.* 108(4):848-858.
5. Holmich, L.R., et al. 2001. Prevalence of silicone breast implant rupture among Danish women. *Plast. Reconstr. Surg.* 108(4):848-858.
6. Collis, N. and D.T. Sharpe. 2000. Silicone gel-filled breast implant integrity: A retrospective review of 478 consecutively explanted implants. *Plast. Reconstr. Surg.* 105:1979-1985.
7. Holmich, L.R., et al. 2004. Untreated silicone breast implant rupture. *Plast. Reconstr. Surg.* 114:204-214.
8. Holmich, L.R., et al. 2004. Untreated silicone breast implant rupture. *Plast. Reconstr. Surg.* 114:204-214.
9. Katzin, W.E., et al. 2005. Pathology of lymph nodes from patients with breast implants: a histologic and spectroscopic evaluation. *Am J Surg Pathol.* 29(4):506-511.
10. Berner, I., M., et al. 2002. Comparative examination of complaints of patients with breast-cancer with and without silicone implants. *Eur. J Obstet. Gynecol. Reprod. Biol.* 102:61-66.
11. Brown, S.L., et al. 2001. Silicone gel breast implant rupture, extracapsular silicone, and health status in a population of women. *J. Rheumatol.* 28:996-1003.
12. Hölmich, L.R., et al. 2003b. Self-reported diseases and symptoms by rupture status among unselected Danish women with cosmetic silicone breast implants. *Plast. Reconstr. Surg.* 111:723-732.
13. Wolfe, F. and J. Anderson. 1999. Silicone filled breast implants and the risk of fibromyalgia and rheumatoid arthritis. *J. Rheumatol.* 26:2025-2028.
14. Hölmich, L.R., et al. 2007. Breast implant rupture and connective tissue disease: A review of the literature. *Plast. Reconstr. Surg.* 120 (Suppl. 1):62S-69S.
15. Hölmich, L.R., et al. 2003b. Self-reported diseases and symptoms by rupture status among unselected Danish women with cosmetic silicone breast implants. *Plast. Reconstr. Surg.* 111:723-732.
16. Bondurant, S., V.L. Ernster and R. Herdman, Eds. 2000. Safety of silicone breast implants. Committee on the Safety of Silicone Breast Implants, Division of Health Promotion and Disease Prevention, Institute of Medicine. Washington, D.C.: National Academy Press.
17. Katzin, W.E., et al. 2005. Pathology of lymph nodes from patients with breast implants: a histologic and spectroscopic evaluation. *Am. J. Surg. Pathol.* 29(4):506-511.
18. Bondurant, S., V.L. Ernster and R. Herdman, Eds. 2000. Safety of silicone breast implants. Committee on the Safety of Silicone Breast Implants, Division of Health Promotion and Disease Prevention, Institute of Medicine. Washington, D.C.: National Academy Press.
19. Brinton, L.A., et al. 2004. Risk of connective tissue disorders among breast implant patients. *Am. J. Epidemiol.* 160(7):619-627.
20. Brown, S.L., et al. 2001. Silicone gel breast implant rupture, extracapsular silicone, and health status in a population of women. *J. Rheumatol.* 28:996-1003.

21. Hölmich, L.R., et al. 2003b. Self-reported diseases and symptoms by rupture status among unselected Danish women with cosmetic silicone breast implants. *Plast. Reconstr. Surg.* 111:723-732.
22. Janowsky, E.C., et al. 2000. Meta-Analyses of the Relation Between Silicone Breast Implants and the Risk of Connective-Tissue Diseases. *N. Engl. J. Med.* 342(11):781-790.
23. Lipworth, L.R.E., et al. 2004. Silicone breast implants and connective tissue disease: An updated review of the epidemiologic evidence. *Ann. Plast. Surg.* 52:598-601.
24. Tugwell, P., et al. 2001. Do silicone breast implants cause rheumatologic disorders? A systematic review for a court-appointed national science panel. *Arthritis Rheum.* (11):2477-2484.
25. Weisman, M.H., et al. 1988. Connective-tissue disease following breast augmentation: A preliminary test of the human adjuvant tissue hypothesis. *Plast. Reconstr. Surg.* 82(4):626-630.
26. Williams, H.J., et al. 1997. Breast implants in patients with differentiated and undifferentiated connective tissue disease. *Arthritis and Rheumatism* 40(3):437-440.
27. Wolfe, F. and J. Anderson. 1999. Silicone filled breast implants and the risk of fibromyalgia and rheumatoid arthritis. *J. Rheumatol.* 26:2025-2028.
28. Bondurant, S., V.L. Ernster and R. Herdman, Eds. 2000. Safety of silicone breast implants. Committee on the Safety of Silicone Breast Implants, Division of Health Promotion and Disease Prevention, Institute of Medicine. Washington, D.C.: National Academy Press.
29. Janowsky, E.C., et al. 2000. Meta-Analyses of the Relation Between Silicone Breast Implants and the Risk of Connective-Tissue Diseases. *N. Engl. J. Med.* 342(11):781-790.
30. Lipworth, L.R.E., et al. 2004. Silicone breast implants and connective tissue disease: An updated review of the epidemiologic evidence. *Ann. Plast. Surg.* 52:598-601.
31. Tugwell, P., et al. 2001. Do silicone breast implants cause rheumatologic disorders? A systematic review for a court-appointed national science panel. *Arthritis Rheum.* (11):2477-2484.
32. Hölmich, L.R., et al. 2003b. Self-reported diseases and symptoms by rupture status among unselected Danish women with cosmetic silicone breast implants. *Plast. Reconstr. Surg.* 111:723-732.
33. Berner, I., M., et al. 2002. Comparative examination of complaints of patients with breast-cancer with and without silicone implants. *Eur. J Obstet. Gynecol. Reprod. Biol.* 102:61-66.
34. Bondurant, S., V.L. Ernster and R. Herdman, Eds. 2000. Safety of silicone breast implants. Committee on the Safety of Silicone Breast Implants, Division of Health Promotion and Disease Prevention, Institute of Medicine. Washington, D.C.: National Academy Press.
35. Breiting, V.B., et al. 2004. Long-term health status of Danish women with silicone breast implants. *Plast. Reconstr. Surg.* 114:217-226.
36. Fryzek, J.P., et al. 2001. Self-reported symptoms among women after cosmetic breast implant and breast reduction surgery. *Plast. Reconstr. Surg.* 107:206-213.
37. Kjølner, K., et al. 2004. Self-reported musculoskeletal symptoms among Danish women with cosmetic breast implants. *Ann Plast Surg.* 52(1):1-7.
38. Brinton, L.A., et al. 2000. Breast cancer following augmentation mammoplasty (United States). *Cancer Causes Control.* 11(9):819-27. *J. Long Term Eff. Med. Implants.* 12(4):271-279.
39. Bryant, H., and Brasher, P. 1995. Breast implants and breast cancer--reanalysis of a linkage study. *N. Engl. J. Med.* 332(23):1535-1539.

40. Deapen, D.M., et al. 1997. Are breast implants anticarcinogenic? A 14-year follow-up of the Los Angeles Study. *Plast. Reconstr. Surg.* 1997 99(5):1346-1353.
41. Herdman, R.C., et al. 2001. Silicone breast implants and cancer. *Cancer Invest.* 2001;19(8):821-832.
42. Pukkala, E., et al. 2002. Incidence of breast and other cancers among Finnish women with cosmetic breast implants, 1970-1999. *J. Long Term Eff. Med. Implants* 12(4):271-279.
43. Brinton, L.A., et al. 2000. Breast cancer following augmentation mammoplasty (United States). *Cancer Causes Control.* 11(9):819-27. *J. Long Term Eff. Med. Implants.* 12(4):271-279.
44. Deapen, D., et al. 2000. Breast cancer stage at diagnosis and survival among patients with prior breast implants. *Plast Reconstr Surg.* 105(2):535-540.
45. Jakubietz, M.G., et al. 2004. Breast augmentation: Cancer concerns and mammography – A literature review. *Plast. Reconstr. Surg.* 113:117e-122e.
46. Miglioretti, D.L., et al. 2004. Effect of breast augmentation on the accuracy of mammography and cancer characteristics. *JAMA* 291(4):442-450.
47. Pukkala, E., et al. 2002. Incidence of breast and other cancers among Finnish women with cosmetic breast implants, 1970-1999. *J. Long Term Eff. Med. Implants* 12(4):271-279.
48. Handel, N. 2007. the effect of silicone implants on the diagnosis, prognosis, and treatment of breast cancer. *Plast. Reconstr. Surg.* 120 (Suppl. 1):81S-93S.
49. Deapen, D. 2007. Breast implants and breast cancer: A review of incidence, detection, mortality, and survival. *Plast. Reconstr. Surg.* 120 (Suppl. 1):70S-80S.
50. Brinton, L.A., et al. 2001b. Cancer risk at sites other than the breast following augmentation mammoplasty. *Ann. Epidemiol.* 11:248-256.
51. McLaughlin, J.K. and L. Lipworth. 2004. Brain cancer and cosmetic breast implants: A review of the epidemiological evidence. *Ann. Plast. Surg.* 52(2):15-17.
52. Brinton, L.A. 2007. The relationship of silicone breast implants and cancer at other sites. *Plast. Reconstr. Surg.* 120 (Suppl. 1):94S-102S.
53. Brinton, L.A., et al. 2001b. Cancer risk at sites other than the breast following augmentation mammoplasty. *Ann. Epidemiol.* 11:248-256.
54. Cook, L.S. 1997. Characteristics of women with and without breast augmentation. *J. Amer. Med. Assoc.* 20:1612-1617.
55. Fryzek, J.P., et al. 2000. Characteristics of women with cosmetic breast augmentation surgery compared with breast reduction surgery patients and women in the general population of Sweden. *Ann Plast Surg.* 45(4):349-356.
56. Kjølner K., et al. 2003. Characteristics of women with cosmetic breast implants compared with women with other types of cosmetic surgery and population-based controls in Denmark. *Ann Plast Surg.* 50(1):6-12.
57. Brinton, L.A., et al. 2001b. Cancer risk at sites other than the breast following augmentation mammoplasty. *Ann. Epidemiol.* 11:248-256.
58. Brinton, L.A. 2007. The relationship of silicone breast implants and cancer at other sites. *Plast. Reconstr. Surg.* 120 (Suppl. 1):94S-102S.
59. Loch-Wilkinson A, Beath KJ, Knight RJW, Wessels WLF, Magnusson M, Papadopoulos T, Connell T, Lofts J, Locke M, Hopper I, Cooter R, Vickery K, Joshi PA, Prince HM, Deva AK. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plast Reconstr Surg.* 2017;140(4):645-654

60. Brody, G.S., et al., Anaplastic Large Cell Lymphoma Occurring in Women with Breast Implants: Analysis of 173 Cases. *Plastic and Reconstructive Surgery*, 2015.135(3):p.695-705.
61. Doren, E.L., et al., U.S. Epidemiology of Breast Implant–Associated Anaplastic Large Cell Lymphoma. *Plastic and Reconstructive Surgery*, 2017. 139(5): p. 1042-1050
62. Brinton, L.A., et al. 2001b. Cancer risk at sites other than the breast following augmentation mammoplasty. *Ann. Epidemiol.* 11:248-256.
63. Bondurant, S., V.L. Ernster and R. Herdman, Eds. 2000. Safety of silicone breast implants. Committee on the Safety of Silicone Breast Implants, Division of Health Promotion and Disease Prevention, Institute of Medicine. Washington, D.C.: National Academy Press.
64. Brinton, L.A., et al. 2001a. Mortality among augmentation mammoplasty patients. *Epidemiol.* 12(3):321-326.
65. Jacobsen, P.H., et al. 2004. Mortality and suicide among Danish women with cosmetic breast implants. *Arch. Int. med.* 164(22):2450-2455.
66. Koot, V., et al. 2003. Total and cost specific mortality among Swedish women with cosmetic breast implants: prospective study. *Br. J. Med.* 326(7388):527-528.
67. Pukkala, E., et al. 2003. Causes of death among Finnish women with cosmetic breast implants, 1971-2001. *Ann. Plast. Surg.* 51(4):339-342.
68. Jacobsen, P.H., et al. 2004. Mortality and suicide among Danish women with cosmetic breast implants. *Arch. Int. med.* 164(22):2450-2455.
69. Lugowski, S.J., et al. 2000. Analysis of silicon in human tissues with special reference to silicone breast implants. *J. Trace Elem. Med. Biol.* 14(1):31-42.
70. Kjølner, K., et al. 2002. Health outcomes in offspring of Danish mothers with cosmetic breast implants. *Ann. Plast. Surg.* 48:238-245.
71. Signorello, L.B., et al. 2001. Offspring health risk after cosmetic breast implantation in Sweden. *Ann. Plast. Surg.* 46:279-286.
72. Hemminki, E., et al. 2004. Births and perinatal health of infants among women who have had silicone breast implantation in Finland, 1967-2000. *Acta Obstet Gynecol Scand.* 83(12):1135-1140.
73. Kjølner, K., et al. 2007. Adverse health outcomes in offspring of mothers with cosmetic breast implants: A review. *Plast. Reconstr. Surg.* 120 (Suppl. 1): 129S-134S.
74. Flassbeck, D.B., et al. 2003. Determination of siloxanes, silicon, and platinum in tissues of women with silicone gel-filled implants. 375(3):356-62 (for example, data from Patients B & C).
75. Bondurant, S., V.L. Ernster and R. Herdman, Eds. 2000. Safety of silicone breast implants. Committee on the Safety of Silicone Breast Implants, Division of Health Promotion and Disease Prevention, Institute of Medicine. Washington, D.C.: National Academy Press.
76. Bondurant, S., V.L. Ernster and R. Herdman, Eds. 2000. Safety of silicone breast implants. Committee on the Safety of Silicone Breast Implants, Division of Health Promotion and Disease Prevention, Institute of Medicine. Washington, D.C.: National Academy Press.
77. Katzin, W.E., et al. 2005. Pathology of lymph nodes from patients with breast implants: a histologic and spectroscopic evaluation. *Am J Surg Pathol.*29(4):506-511.

- 
78. Stein, J., et al. 1999. In situ determination of the active catalyst in hydrosilylation reactions using highly reactive Pt(0) catalyst precursors. *J. Am. Chem. Soc.* 121(15):3693-3703. Chandra, G., et al. 1987. A convenient and novel route to bis(alkyne)platinum(0) and other platinum(0) complexes from Speier's hydrosilylation catalyst. *Organometallics*. 6:191-2. Lappert, M.F. and Scott, F.P.A. 1995. The reaction pathway from Speier's to Karstedt's hydrosilylation catalyst. *J. Organomet. Chem.* 492(2):C11-C13. Lewis, L.N., et al. 1995. Mechanism of formation of platinum(0) complexes containing silicon-vinyl ligands. *Organometallics*. 14:2202-2213.
  79. Henriksen, T.F., et al. 2005. Surgical intervention and capsular contracture after breast augmentation: a prospective study of risk factors. *Ann. Plast. Surg.* 54(4):343-351.
  80. For example: Seify, H., et al. 2005. Preliminary (3 years) experience with smooth wall silicone gel implants for primary breast augmentation. 54(3):231-235.
  81. Henriksen, T.F., et al. 2005. Surgical intervention and capsular contracture after breast augmentation: a prospective study of risk factors. *Ann. Plast. Surg.* 54(4):343-351.
  82. For example: Kulmala, I., et al. 2004. Local complications after cosmetic breast implant surgery in Finland. *Ann. Plast. Surg.* 53(5):413-419.
  83. Bondurant, S., V.L. Ernster and R. Herdman, Eds. 2000. Safety of silicone breast implants. Committee on the Safety of Silicone Breast Implants, Division of Health Promotion and Disease Prevention, Institute of Medicine Washington, D.C.: National Academy Press.
  84. Henriksen, T.F., et al. 2005. Surgical intervention and capsular contracture after breast augmentation: a prospective study of risk factors. *Ann. Plast. Surg.* 54(4):343-351.
  85. For example: Kulmala, I., et al. 2004. Local complications after cosmetic breast implant surgery in Finland. *Ann. Plast. Surg.* 53(5):413-419.
  86. Hurst, N.M. 1996. Lactation after augmentation mammoplasty. *Obstet. Gynecol.* 87(1):30-34.

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40. Deapen, D.M., et al. 1997. Are breast implants anticarcinogenic? A 14-year follow-up of the Los Angeles Study. *Plast. Reconstr. Surg.* 1997 99(5):1346-1353.
41. Herdman, R.C., et al. 2001. Silicone breast implants and cancer. *Cancer Invest.* 2001;19(8):821-832.
42. Pukkala, E., et al. 2002. Incidence of breast and other cancers among Finnish women with cosmetic breast implants, 1970-1999. *J. Long Term Eff. Med. Implants* 12(4):271-279.
43. Brinton, L.A., et al. 2000. Breast cancer following augmentation mammoplasty (United States). *Cancer Causes Control.* 11(9):819-27. *J. Long Term Eff. Med. Implants.* 12(4):271-279.
44. Deapen, D., et al. 2000. Breast cancer stage at diagnosis and survival among patients with prior breast implants. *Plast Reconstr Surg.* 105(2):535-540.
45. Jakubietz, M.G., et al. 2004. Breast augmentation: Cancer concerns and mammography – A literature review. *Plast. Reconstr. Surg.* 113:117e-122e.
46. Miglioretti, D.L., et al. 2004. Effect of breast augmentation on the accuracy of mammography and cancer characteristics. *JAMA* 291(4):442-450.
47. Pukkala, E., et al. 2002. Incidence of breast and other cancers among Finnish women with cosmetic breast implants, 1970-1999. *J. Long Term Eff. Med. Implants* 12(4):271-279.
48. Handel, N. 2007. the effect of silicone implants on the diagnosis, prognosis, and treatment of breast cancer. *Plast. Reconstr. Surg.* 120 (Suppl. 1):81S-93S.
49. Deapen, D. 2007. Breast implants and breast cancer: A review of incidence, detection, mortality, and survival. *Plast. Reconstr. Surg.* 120 (Suppl. 1):70S-80S.
50. Brinton, L.A., et al. 2001b. Cancer risk at sites other than the breast following augmentation mammoplasty. *Ann. Epidemiol.* 11:248-256.
51. McLaughlin, J.K. and L. Lipworth. 2004. Brain cancer and cosmetic breast implants: A review of the epidemiological evidence. *Ann. Plast. Surg.* 52(2):15-17.
52. Brinton, L.A. 2007. The relationship of silicone breast implants and cancer at other sites. *Plast. Reconstr. Surg.* 120 (Suppl. 1):94S-102S.
53. Brinton, L.A., et al. 2001b. Cancer risk at sites other than the breast following augmentation mammoplasty. *Ann. Epidemiol.* 11:248-256.
54. Cook, L.S. 1997. Characteristics of women with and without breast augmentation. *J. Amer. Med. Assoc.* 20:1612-1617.
55. Fryzek, J.P., et al. 2000. Characteristics of women with cosmetic breast augmentation surgery compared with breast reduction surgery patients and women in the general population of Sweden. *Ann Plast Surg.* 45(4):349-356.
56. Kjølner K., et al. 2003. Characteristics of women with cosmetic breast implants compared with women with other types of cosmetic surgery and population-based controls in Denmark. *Ann Plast Surg.* 50(1):6-12.
57. Brinton, L.A., et al. 2001b. Cancer risk at sites other than the breast following augmentation mammoplasty. *Ann. Epidemiol.* 11:248-256.
58. Brinton, L.A. 2007. The relationship of silicone breast implants and cancer at other sites. *Plast. Reconstr. Surg.* 120 (Suppl. 1):94S-102S.
59. Loch-Wilkinson A, Beath KJ, Knight RJW, Wessels WLF, Magnusson M, Papadopoulos T, Connell T, Lofts J, Locke M, Hopper I, Cooter R, Vickery K, Joshi PA, Prince HM, Deva AK. Breast Implant-Associated Anaplastic Large Cell Lymphoma in Australia and New Zealand: High-Surface-Area Textured Implants Are Associated with Increased Risk. *Plast Reconstr Surg.* 2017;140(4):645-654.

60. Brody, G.S., et al., Anaplastic Large Cell Lymphoma Occurring in Women with Breast Implants: Analysis of 173 Cases. *Plastic and Reconstructive Surgery*, 2015.135(3):p.695-705.
61. Doren, E.L., et al., U.S. Epidemiology of Breast Implant–Associated Anaplastic Large Cell Lymphoma. *Plastic and Reconstructive Surgery*, 2017. 139(5): p. 1042-1050.
62. Brinton, L.A., et al. 2001b. Cancer risk at sites other than the breast following augmentation mammoplasty. *Ann. Epidemiol.* 11:248-256.
63. Bondurant, S., V.L. Ernster and R. Herdman, Eds. 2000. Safety of silicone breast implants. Committee on the Safety of Silicone Breast Implants, Division of Health Promotion and Disease Prevention, Institute of Medicine. Washington, D.C.: National Academy Press.
64. Brinton, L.A., et al. 2001a. Mortality among augmentation mammoplasty patients. *Epidemiol.* 12(3):321-326.
65. Jacobsen, P.H., et al. 2004. Mortality and suicide among Danish women with cosmetic breast implants. *Arch. Int. med.* 164(22):2450-2455.
66. Koot, V., et al. 2003. Total and cost specific mortality among Swedish women with cosmetic breast implants: prospective study. *Br. J. Med.* 326(7388):527-528.
67. Pukkala, E., et al. 2003. Causes of death among Finnish women with cosmetic breast implants, 1971-2001. *Ann. Plast. Surg.* 51(4):339-342.
68. Jacobsen, P.H., et al. 2004. Mortality and suicide among Danish women with cosmetic breast implants. *Arch. Int. med.* 164(22):2450-2455.
69. Lugowski, S.J., et al. 2000. Analysis of silicon in human tissues with special reference to silicone breast implants. *J. Trace Elem. Med. Biol.* 14(1):31-42.
70. Kjølner, K., et al. 2002. Health outcomes in offspring of Danish mothers with cosmetic breast implants. *Ann. Plast. Surg.* 48:238-245.
71. Signorello, L.B., et al. 2001. Offspring health risk after cosmetic breast implantation in Sweden. *Ann. Plast. Surg.* 46:279-286.
72. Hemminki, E., et al. 2004. Births and perinatal health of infants among women who have had silicone breast implantation in Finland, 1967-2000. *Acta Obstet Gynecol Scand.* 83(12):1135-1140.
73. Kjølner, K., et al. 2007. Adverse health outcomes in offspring of mothers with cosmetic breast implants: A review. *Plast. Reconstr. Surg.* 120 (Suppl. 1): 129S-134S.
74. Flassbeck, D.B., et al. 2003. Determination of siloxanes, silicon, and platinum in tissues of women with silicone gel-filled implants. 375(3):356-62 (for example, data from Patients B & C).
75. Bondurant, S., V.L. Ernster and R. Herdman, Eds. 2000. Safety of silicone breast implants. Committee on the Safety of Silicone Breast Implants, Division of Health Promotion and Disease Prevention, Institute of Medicine. Washington, D.C.: National Academy Press.
76. Bondurant, S., V.L. Ernster and R. Herdman, Eds. 2000. Safety of silicone breast implants. Committee on the Safety of Silicone Breast Implants, Division of Health Promotion and Disease Prevention, Institute of Medicine. Washington, D.C.: National Academy Press.
77. Katzin, W.E., et al. 2005. Pathology of lymph nodes from patients with breast implants: a histologic and spectroscopic evaluation. *Am J Surg Pathol.*29(4):506-511.

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**For customer service call** Mentor, a unit of Johnson & Johnson Medical Products, a division of Johnson & Johnson Inc., **at** 1-800-668-6069 **or contact your local Mentor representative.**  
[www.mentorwllc.com](http://www.mentorwllc.com) • [www.mentordirect.com](http://www.mentordirect.com)

**Manufacturer**

MENTOR  
3041 Skyway Circle North  
Irving, TX 75038-3540  
USA  
972-252-6060

**Manufacturer**

Mentor Medical Systems B.V.  
Zernikedreef 2,  
2333 CL Leiden  
Netherlands  
+31 71 751360

\* For recognized manufacturer, refer to product label