

LCP Proximal Radius Plates 2.4

Plates for internal fixation of the proximal radius fractures

Surgical Technique



 Image intensifier control

This description alone does not provide sufficient background for direct use of DePuy Synthes products. Instruction by a surgeon experienced in handling these products is highly recommended.

Processing, Reprocessing, Care and Maintenance

For general guidelines, function control and dismantling of multi-part instruments, as well as processing guidelines for implants, please contact your local sales representative or refer to:

<http://emea.depuySynthes.com/hcp/reprocessing-care-maintenance>

For general information about reprocessing, care and maintenance of DePuy Synthes reusable devices, instrument trays and cases, as well as processing of DePuy Synthes non-sterile implants, please consult the Important Information leaflet (SE_023827) or refer to:

<http://emea.depuySynthes.com/hcp/reprocessing-care-maintenance>

Table of Contents

| | | |
|----------------------------|---|----|
| Introduction | LCP Proximal Radius Plates 2.4 | 2 |
| | The AO Principles of Fracture Management | 3 |
| | | |
| Surgical Technique | Preparation | 4 |
| | Reduction | 6 |
| | Proximal Screw Insertion | 7 |
| | Shaft Screw Insertion | 8 |
| | General Notes on Technique | 9 |
| | Postoperative Treatment and Implant Removal | 12 |
| | | |
| Product Information | Plates | 13 |
| | Screws | 14 |
| | Instruments | 15 |
| | | |
| MRI Information | | 16 |

- **Notes**
- ▲ **Precautions**
- ▲ **WARNINGS**

LCP Proximal Radius Plates 2.4

Plate features

- Nine LCP Proximal Radius Plates available to address internal fixation of the proximal radius fractures
- Plates are contoured
- Combi holes allow fixation with locking screws in the threaded section and cortex screws in the Dynamic Compression Unit (DCU) section for distraction.

Additional features

- Limited-contact design shaft with 2, 3, and 4 combi-holes
- The holes in the head of the plate accept 2.4 mm locking screws
- The shaft holes accept 2.4 mm locking screws in the threaded portion or 2.7 mm cortex screws and 2.4 mm cortex screws in the distraction portion
- Plates for radial head rim available in right and left plates with a 5° tilt
- Plates for radial head neck fit both the left and right side of the proximal radius



Intended Use, Indications and Contraindications can be found in the corresponding system Instructions for Use.

The AO Principles of Fracture Management

Mission


The AO's mission is promoting excellence in patient care and outcomes in trauma and musculoskeletal disorders.

AO Principles^{1,2}

1.




2.



3.



4.



Fracture reduction and fixation to restore anatomical relationships.

Fracture fixation providing absolute or relative stability, as required by the "personality" of the fracture, the patient, and the injury.

Preservation of the blood supply to soft-tissues and bone by gentle reduction techniques and careful handling.

Early and safe mobilization and rehabilitation of the injured part and the patient as a whole.


¹ Müller ME, Allgöwer M, Schneider R, Willenegger H. Manual of Internal Fixation. 3rd ed. Berlin, Heidelberg New York: Springer 1991.

² Buckley RE, Moran CG, Apivatthakakul T. AO Principles of Fracture Management: 3rd ed. Vol. 1: Principles, Vol. 2: Specific fractures. Thieme; 2017.

Preparation


Patient position

Position the patient supine on the table. Prepare the extremity from the axilla to the hand. This allows rotation of the forearm and flexion and extension of the elbow during the operative fixation.



Approach


A lateral approach is most commonly used. Take care to avoid the deep branch of the radial nerve, which runs anterior to the capsule and the radial head. To minimize the risk of operative disruption of the lateral collateral ligament, the capsular incision should remain in the front of the anterior margin of the anconeus muscle and parallel to the fascial limit of the extensor carpi ulnaris. The annular ligament, a true thickening of the capsule, is opened laterally or slightly anteriorly to allow full inspection of the fragments. In selected cases, an osteotomy of the lateral epicondyle will allow an extensile approach.



Illustrations on this page are reproduced from Rüedi/Murphy, AO Principles of Fracture Management by permission of AO Publishing, Copyright © 2000 AO Publishing, Switzerland.


Safe zone

Both plates are placed within the Hotchkiss safe zone¹ which is described as an area of 110° on the radial head that is free of impingement between ulna and radius. The Hotchkiss safe zone is located on the opposite side of the radial tuberosity.



Neutral

Supination



Pronation

¹ Smith GR, Hotchkiss RN. Radial head and neck fractures: anatomic guidelines for proper placement of internal fixation. J Shoulder Elbow Surg. 1996 Mar-Apr;5(2 Pt 1):113-7.

Reduction

1. Reduce fracture

Perform preliminary fixation with Kirschner wires.

▲ **Precaution:**

The plate reduction wires and Kirschner wires are single use items, do not re-use.



2. Apply Radial Head Plate

Instruments

311.430 Handle with Quick Coupling,
length 110 mm

314.467 Screwdriver Shaft, Stardrive, T8,
self-holding

After reducing the fracture, apply the plate and insert a preliminary cortex screw through the DCU portion of one of the elongated Combi holes.

■ **Note:**

The plates are contoured. There are no undercuts on the plates to protect the threaded holes from distortion during bending. If it is necessary to bend the plate, use bending pliers.

▲ **Precaution:**

Reverse bending or use of the incorrect instrumentation for bending may weaken the plate and lead to premature plate failure (e.g. breakage). Do not bend the plate beyond what is required to match the anatomy.



Proximal Screw Insertion

3. Insert proximal screws

Instruments

| | |
|------------|---|
| 311.430 | Handle with Quick Coupling, length 110 mm |
| 03.110.005 | Handle for Torque Limiters 0.4/0.8/1.2 Nm |
| 314.467 | Screwdriver Shaft, Stardrive, T8, self-holding |
| 511.776 | Torque Limiter 0.8 Nm, with Quick Coupling |
| 310.509 | Drill Bit Ø 1.8 mm, with marking, length 110/85 mm, 2-flute, for Quick Coupling |
| 323.029 | LCP Drill Sleeve 2.4, with Scale up to 30 mm, for Drill Bits Ø 1.8 mm |

Use 2.4 mm locking screws in the head of the plate.

▲ Precaution:

Careful drilling is necessary, as interference with screws in the proximal portion of the plate is possible. In case of interference, stop drilling and use a screw of appropriate length.

See General Notes on Technique for additional information.



Shaft Screw Insertion

4. Insert screws in the shaft of the plate

Instruments

311.430 Handle with Quick Coupling,
length 110 mm

314.467 Screwdriver Shaft, Stardrive, T8,
self-holding

For cortex screws 2.4 mm

310.509 Drill Bit Ø 1.8 mm, with marking,
length 110/85 mm, 2-flute,
for Quick Coupling

310.530 Drill Bit Ø 2.4 mm, length 100/75 mm,
2-flute, for Quick Coupling

323.202 Universal Drill Guide 2.4

For cortex screws 2.7 mm

310.534 Drill Bit Ø 2.0 mm, length 100/75 mm,
2-flute, for Quick Coupling

310.260 Drill Bit Ø 2.7 mm, length 100/75 mm,
2-flute, for Quick Coupling

323.260 Universal Drill Guide 2.7

For locking screws 2.4 mm

310.509 Drill Bit Ø 1.8 mm, with marking,
length 110/85 mm, 2-flute,
for Quick Coupling

323.029 LCP Drill Sleeve 2.4, with Scale up
to 30 mm, for Drill Bits Ø 1.8 mm

Use 2.4 mm and 2.7 mm cortex screws and 2.4 mm
locking screws in the shaft of the plate.



General Notes on Technique

Determine screw choice

If planning a combination of locking and cortex screws, first use a cortex screw to pull the plate to the bone.

▲ Precaution:

If using a locking screw first, ensure that the plate is held securely to the bone, to avoid spinning of the plate as the screw locks into the plate.

Insert cortex screws

Instruments

| | |
|---------|---|
| 310.534 | Drill Bit Ø 2.0 mm, with marking, length 110/85 mm, 2-flute, for Quick Coupling |
| 310.260 | Drill Bit Ø 2.7 mm, length 100/75 mm, 2-flute, for Quick Coupling |
| 310.509 | Drill Bit Ø 1.8 mm, with marking, length 110/85 mm, 2-flute, for Quick Coupling |
| 310.530 | Drill Bit Ø 2.4 mm, length 100/75 mm, 2-flute, for Quick Coupling |
| 319.005 | Depth Gauge for Screws Ø 2.0 and 2.4 mm, measuring range up to 40 mm |
| 323.202 | Universal Drill Guide 2.4 |
| 323.260 | Universal Drill Guide 2.7 |

- Use the 2.4 or 2.7 universal drill guide for an eccentric (distraction) or neutral (butress) insertion of cortex screws.
- For the 2.4 mm cortex screw, use the 1.8 mm drill bit for the threaded hole and the 2.4 mm drill bit for the gliding hole. Determine the length of the screw by using the depth gauge. For 2.7 mm cortex screws, use the 2.0 mm drill bit for the threaded hole and the 2.7 mm drill bit for the gliding hole.

Insert locking screws

Instruments

| | |
|------------|---|
| 310.509 | Drill Bit Ø 1.8 mm, with marking, length 110/85 mm, 2-flute, for Quick Coupling |
| 311.430 | Handle with Quick Coupling, length 110 mm |
| 03.110.005 | Handle for Torque Limiters 0.4/0.8/1.2 Nm |
| 314.467 | Screwdriver Shaft, Stardrive, T8, self-holding |
| 511.776 | Torque Limiter 0.8 Nm, with Quick Coupling |
| 319.005 | Depth Gauge for Screws Ø 2.0 and 2.4 mm, measuring range up to 40 mm |
| 323.029 | LCP Drill Sleeve 2.4, with Scale up to 30 mm, for Drill Bits Ø 1.8 mm |



- Screw the LCP drill guide into a hole until fully seated.
- For the 2.4 mm locking screws, use the 1.8 mm drill bit with marking, to drill to the desired depth. Determine the screw length directly from the mark on the drill bit and scale on the drill sleeve. This may be verified using a depth gauge for 2.0 mm and 2.4 mm screws.
- Insert the locking screw manually with the screwdriver shaft handle and Torque Limiter 0.8 Nm. Carefully tighten the locking screw.

Alternative: Locking screw insertion with holding sleeve

Instruments

| | |
|------------|--|
| 311.430 | Handle with Quick Coupling, length 110 mm |
| 03.110.005 | Handle for Torque Limiters 0.4/0.8/1.2 Nm |
| 314.467 | Screwdriver Shaft, Stardrive, T8, self-holding |
| 511.776 | Torque Limiter 0.8 Nm, with Quick Coupling |
| 314.468 | Holding Sleeve for Screws Stardrive Ø 2.4 mm, T8, for Screwdriver Shafts Ø 3.5 mm, for 314.467 |



An alternate method may be used for insertion of locking screws, using the locking screw to pull the plate to the bone. Place the holding sleeve onto the screwdriver shaft. Pick up the locking screw with the holding sleeve and insert the screw. With the locking screw still held by the holding sleeve, tighten the screw until the plate is drawn to the bone. Pull up on the holding sleeve to release the screw head, and tighten the screw until a click is heard.

▲ Precaution:

Take care to hold the plate securely on the bone to avoid spinning the plate as the screw locks to the plate.

Postoperative Treatment and Implant Removal

Postoperative treatment

Postoperative treatment with locking compression plates does not differ from conventional internal fixation procedures.

Implant removal

To remove locking screws, first unlock all screws from the plate and then remove the screws completely from the bone. This avoids rotation of the plate when removing the last locking screw.

Plates

Rim plates

2-hole plates X41.680 (right)

X41.681 (left)

3-hole plates X41.682 (right)

X41.683 (left)

4-hole plates X41.684 (right)

X41.685 (left)



Neck plates

2-hole plates X41.690

3-hole plates X41.691

4-hole plates X41.692



All implants are available non-sterile or sterile packed.
Add suffix "S" to article number to order sterile product.

X = 2: Stainless Steel
X = 4: Titan (TiCP)

Screws

Cortex Screws 2.4 mm, self-tapping, Stardrive, T8

(X01.760 – X01.780, length 10–30 mm)

- For use in Combi holes in plate shaft
- Low-profile head in the plate holes
- Stardrive recess mates with self-retaining screwdriver
- 10 mm to 30 mm lengths



Locking Screws 2.4 mm, self-tapping, Stardrive, T8

(X12.806 – X12.830, length 6–30 mm)

- Threaded, conical head locks securely into the plate
- Locked screws facilitate unicortical screw fixation
- Stardrive recess mates with self-retaining screwdriver
- 6 mm to 30 mm lengths



Cortex Screws 2.7 mm, self-tapping, Stardrive, T8

(X02.870 – X02.886, length 10–26 mm)

- For use in Combi holes in plate shaft
- Low-profile head in the plate holes
- Stardrive recess mates with self-retaining screwdriver
- 10 mm to 26 mm lengths




■ Note:

For information on fixation principles using conventional and locked plating techniques, please refer to the LCP Surgical Technique.

All implants are available non-sterile or sterile packed.
Add suffix "S" to article number to order sterile product.

X = 2: Stainless Steel
X = 4: Titanium Alloy (TAN)

Instruments

| | | |
|------------|---|--|
| 311.430 | Handle with Quick Coupling, length 110 mm |  |
| 314.467 | Screwdriver Shaft Stardrive, T8, self-holding |  |
| 314.468 | Holding Sleeve for Screws Stardrive 2.4, for Screwdriver Shaft 314.467 |  |
| 323.029 | LCP Drill Sleeve 2.4, with scale up to 30 mm, for Drill Bit 1.8 mm |  |
| 310.534 | Drill Bit Ø 2.0 mm with marking, length 110/85 mm, 2-fluted, for Quick Coupling |  |
| 310.509 | Drill Bit Ø 1.8 mm with marking, length 110/85 mm, 2-fluted, for Quick Coupling |  |
| 511.776 | Torque Limiter 0.8 Nm, with Quick Coupling |  |
| 03.110.005 | Handle for Torque Limiters 0.4/0.8/1.2 Nm |  |

MRI Information

Torque, Displacement and Image Artifacts according to ASTM F 2213, ASTM F 2052 and ASTM F 2119

Non-clinical testing of worst case scenario in a 3 T MRI system did not reveal any relevant torque or displacement of the construct for an experimentally measured local spatial gradient of the magnetic field of 3.69 T/m. The largest image artifact extended approximately 169 mm from the construct when scanned using the Gradient Echo (GE). Testing was conducted on a 3 T MRI system.

Radio-Frequency-(RF-)induced heating according to ASTM F 2182

Non-clinical electromagnetic and thermal simulations of a worst case scenario lead to temperature rises of 7.5 °C (1.5 T) and 6.1 °C (3 T) under MRI Conditions using RF Coils (whole body averaged specific absorption rate [SAR] of 2 W/kg for 15 minutes).

▲ Precautions:

The above mentioned test relies on non-clinical testing. The actual temperature rise in the patient will depend on a variety of factors beyond the SAR and time of RF application. Thus, it is recommended to pay particular attention to the following points:

- It is recommended to thoroughly monitor patients undergoing MR scanning for perceived temperature and/or pain sensations.
- Patients with impaired thermoregulation or temperature sensation should be excluded from MR scanning procedures.
- Generally, it is recommended to use a MR system with low field strength in the presence of conductive implants. The employed specific absorption rate (SAR) should be reduced as far as possible.
- Using the ventilation system may further contribute to reduce temperature increase in the body.

Not all products are currently available in all markets.
This publication is not intended for distribution in the USA.
Intended use, Indications and Contraindications can be found in the corresponding system Instructions for Use.
All Surgical Techniques are available as PDF files at www.depuysynthes.com/ifu



 Synthes GmbH
Eimattstrasse 3
4436 Oberdorf
Switzerland
Tel: +41 61 965 61 11

www.depuysynthes.com