

Depuy Synthes GmbH
OrthoKit® Rental Service

Im Brunnenfeld 8
D-79224 Umkirch



Tel. +49 800 201 20 70
Fax +49 7665 503 290
www.depuysynthes.com




L004344S

Db. kopfl. KomprSchr. Ø4.5 & 5.5 Steril

Box01/02

Pos.	Mat.Nr	Description	Qty	Unit	Location
Box01/02		Tray 01			
1	03.333.103S	SprBohr Ø 3.0 dbohrt L215 SKpl 	1	ST	B01T01
2	03.333.203S	Kopfraumfräser Ø 4.5/5.5 dbohrt SKpl 	1	ST	B01T01
3	04.333.520S	CCHS - Ø 4.5, L 20mm, Kurzes Gew, TAV	2	ST	B01T01
4	04.333.522S	CCHS - Ø 4.5, L 22mm, Kurzes Gew, TAV	2	ST	B01T01
5	04.333.524S	CCHS - Ø 4.5, L 24mm, Kurzes Gew, TAV	2	ST	B01T01
6	04.333.526S	CCHS - Ø 4.5, L 26mm, Kurzes Gew, TAV	2	ST	B01T01
7	04.333.528S	CCHS - Ø 4.5, L 28mm, Kurzes Gew, TAV	2	ST	B01T01
8	04.333.530S	CCHS - Ø 4.5, L 30mm, Kurzes Gew, TAV	2	ST	B01T01
9	04.333.532S	CCHS - Ø 4.5, L 32mm, Kurzes Gew, TAV	2	ST	B01T01
10	04.333.534S	CCHS - Ø 4.5, L 34mm, Kurzes Gew, TAV	2	ST	B01T01
11	04.333.536S	CCHS - Ø 4.5, L 36mm, Kurzes Gew, TAV	2	ST	B01T01
12	04.333.538S	CCHS - Ø 4.5, L 38mm, Kurzes Gew, TAV	2	ST	B01T01
13	04.333.540S	CCHS - Ø 4.5, L 40mm, Kurzes Gew, TAV	2	ST	B01T01
14	04.333.542S	CCHS - Ø 4.5, L 42mm, Kurzes Gew, TAV	2	ST	B01T01
15	04.333.544S	CCHS - Ø 4.5, L 44mm, Kurzes Gew, TAV	2	ST	B01T01
16	04.333.546S	CCHS - Ø 4.5, L 46mm, Kurzes Gew, TAV	2	ST	B01T01
17	04.333.548S	CCHS - Ø 4.5, L 48mm, Kurzes Gew, TAV	2	ST	B01T01
18	04.333.550S	CCHS - Ø 4.5, L 50mm, Kurzes Gew, TAV	2	ST	B01T01
19	04.333.555S	CCHS - Ø 4.5, L 55mm, Kurzes Gew, TAV	2	ST	B01T01
20	04.333.560S	CCHS - Ø 4.5, L 60mm, Kurzes Gew, TAV	2	ST	B01T01
21	04.333.565S	CCHS - Ø 4.5, L 65mm, Kurzes Gew, TAV	2	ST	B01T01
22	04.333.570S	CCHS - Ø 4.5, L 70mm, Kurzes Gew, TAV	2	ST	B01T01
23	04.333.575S	CCHS - Ø 4.5, L 75mm, Kurzes Gew, TAV	2	ST	B01T01
24	04.333.580S	CCHS - Ø 4.5, L 80mm, Kurzes Gew, TAV	2	ST	B01T01
25	04.333.585S	CCHS - Ø 4.5, L 85mm, Kurzes Gew, TAV	2	ST	B01T01



Depuy Synthes GmbH
OrthoKit® Rental Service

Im Brunnenfeld 8
D-79224 Umkirch

Tel. +49 800 201 20 70
Fax +49 7665 503 290
www.depuyorthosynthes.com



L004344S

Db. kopfl. KomprSchr. Ø4.5 & 5.5 Steril

26	04.333.590S	CCHS - Ø 4.5, L 90mm, Kurzes Gew, TAV		2	ST	B01T01
27	04.333.595S	CCHS - Ø 4.5, L 95mm, Kurzes Gew, TAV		2	ST	B01T01
28	04.333.500S	CCHS - Ø 4.5, L 100mm, Kurzes Gew, TAV		2	ST	B01T01
29	04.333.501S	CCHS - Ø 4.5, L 105mm, Kurzes Gew, TAV		2	ST	B01T01
30	04.333.502S	CCHS - Ø 4.5, L 110mm, Kurzes Gew, TAV		2	ST	B01T01
31	04.333.600S	CHS - Ø 5.5, L 100mm, Kurzes Gew, TAV		2	ST	B01T01
32	04.333.601S	CCHS - Ø 5.5, L 105mm, Kurzes Gew, TAV		2	ST	B01T01
33	04.333.602S	CCHS - Ø 5.5, L 110mm, Kurzes Gew, TAV		2	ST	B01T01
34	04.333.620S	CCHS - Ø 5.5, L 20mm, Kurzes Gew, TAV		2	ST	B01T01
35	04.333.622S	CCHS - Ø 5.5, L 22mm, Kurzes Gew, TAV		2	ST	B01T01
36	04.333.624S	CCHS - Ø 5.5, L 24mm, Kurzes Gew, TAV		2	ST	B01T01
37	04.333.626S	CCHS - Ø 5.5, L 26mm, Kurzes Gew, TAV		2	ST	B01T01
38	04.333.628S	CCHS - Ø 5.5, L 28mm, Kurzes Gew, TAV		2	ST	B01T01
39	04.333.630S	CCHS - Ø 5.5, L 30mm, Kurzes Gew, TAV		2	ST	B01T01
40	04.333.632S	CCHS - Ø 5.5, L 32mm, Kurzes Gew, TAV		2	ST	B01T01
41	04.333.634S	CCHS - Ø 5.5, L 34mm, Kurzes Gew, TAV		2	ST	B01T01
42	04.333.636S	CCHS - Ø 5.5, L 36mm, Kurzes Gew, TAV		2	ST	B01T01
43	04.333.638S	CCHS - Ø 5.5, L 38mm, Kurzes Gew, TAV		2	ST	B01T01
44	04.333.640S	CCHS - Ø 5.5, L 40mm, Kurzes Gew, TAV		2	ST	B01T01
45	04.333.642S	CCHS - Ø 5.5, L 42mm, Kurzes Gew, TAV		2	ST	B01T01
46	04.333.644S	CCHS - Ø 5.5, L 44mm, Kurzes Gew, TAV		2	ST	B01T01
47	04.333.646S	CCHS - Ø 5.5, L 46mm, Kurzes Gew, TAV		2	ST	B01T01
48	04.333.648S	CCHS - Ø 5.5, L 48mm, Kurzes Gew, TAV	2	ST	B01T01	
49	04.333.650S	CCHS - Ø 5.5, L 50mm, Kurzes Gew, TAV	2	ST	B01T01	
50	04.333.655S	CCHS - Ø 5.5, L 55mm, Kurzes Gew, TAV	2	ST	B01T01	
51	04.333.660S	CCHS - Ø 5.5, L 60mm, Kurzes Gew, TAV	2	ST	B01T01	
52	04.333.665S	CCHS - Ø 5.5, L 65mm, Kurzes Gew, TAV	2	ST	B01T01	

Depuy Synthes GmbH Im Brunnenfeld 8
OrthoKit® Rental Service D-79224 Umkirch

Tel. +49 800 201 20 70
Fax +49 7665 503 290
www.depuyorthosynthes.com





L004344S

Db. kopfl. KomprSchr. Ø4.5 & 5.5 Steril

Box02/02

Pos.	Mat.Nr	Description	Qty	Unit	Location
Box02/02		Tray 01			
53	04.333.670S	CCHS - Ø 5.5, L 70mm, Kurzes Gew, TAV	2	ST	B02T01
54	04.333.675S	CCHS - Ø 5.5, L 75mm, Kurzes Gew, TAV	2	ST	B02T01
55	04.333.680S	CCHS - Ø 5.5, L 80mm, Kurzes Gew, TAV	2	ST	B02T01
56	04.333.685S	CCHS - Ø 5.5, L 85mm, Kurzes Gew, TAV	2	ST	B02T01
57	04.333.690S	CCHS - Ø 5.5, L 90mm, Kurzes Gew, TAV	2	ST	B02T01
58	04.333.695S	CCHS - Ø 5.5, L 95mm, Kurzes Gew, TAV	2	ST	B02T01
59	04.334.530S	CCHS - Ø 4.5, L 30mm, Langes Gew, TAV	2	ST	B02T01
60	04.334.532S	CCHS - Ø 4.5, L 32mm, Langes Gew, TAV	2	ST	B02T01
61	04.334.534S	CCHS - Ø 4.5, L 34mm, Langes Gew, TAV	2	ST	B02T01
62	04.334.536S	CCHS - Ø 4.5, L 36mm, Langes Gew, TAV	2	ST	B02T01
63	04.334.538S	CCHS - Ø 4.5, L 38mm, Langes Gew, TAV	2	ST	B02T01
64	04.334.540S	CCHS - Ø 4.5, L 40mm, Langes Gew, TAV	2	ST	B02T01
65	04.334.542S	CCHS - Ø 4.5, L 42mm, Langes Gew, TAV	2	ST	B02T01
66	04.334.544S	CCHS - Ø 4.5, L 44mm, Langes Gew, TAV	2	ST	B02T01
67	04.334.546S	CCHS - Ø 4.5, L 46mm, Langes Gew, TAV	2	ST	B02T01
68	04.334.548S	CCHS - Ø 4.5, L 48mm, Langes Gew, TAV	2	ST	B02T01
69	04.334.550S	CCHS - Ø 4.5, L 50mm, Langes Gew, TAV	2	ST	B02T01
70	04.334.555S	CCHS - Ø 4.5, L 55mm, Langes Gew, TAV	2	ST	B02T01
71	04.334.560S	CCHS - Ø 4.5, L 60mm, Langes Gew, TAV	2	ST	B02T01
72	04.334.565S	CCHS - Ø 4.5, L 65mm, Langes Gew, TAV	2	ST	B02T01
73	04.334.570S	CCHS - Ø 4.5, L 70mm, Langes Gew, TAV	2	ST	B02T01
74	04.334.575S	CCHS - Ø 4.5, L 75mm, Langes Gew, TAV	2	ST	B02T01
75	04.334.580S	CCHS - Ø 4.5, L 80mm, Langes Gew, TAV	2	ST	B02T01
76	04.334.585S	CCHS - Ø 4.5, L 85mm, Langes Gew, TAV	2	ST	B02T01
77	04.334.590S	CCHS - Ø 4.5, L 90mm, Langes Gew, TAV	2	ST	B02T01
78	04.334.595S	CCHS - Ø 4.5, L 95mm, Langes Gew, TAV	2	ST	B02T01



Depuy Synthes GmbH
OrthoKit® Rental Service

Im Brunnenfeld 8
D-79224 Umkirch

Tel. +49 800 201 20 70
Fax +49 7665 503 290
www.depuysynthes.com



L004344S

Db. kopfl. KomprSchr. Ø4.5 & 5.5 Steril

79	04.334.500S	CCHS - Ø 4.5, L 100mm, Langes Gew, TAV		2	ST	B02T01
80	04.334.501S	CCHS - Ø 4.5, L 105mm, Langes Gew, TAV		2	ST	B02T01
81	04.334.502S	CCHS - Ø 4.5, L 110mm, Langes Gew, TAV		2	ST	B02T01
82	04.334.630S	CCHS - Ø 5.5, L 30mm, Langes Gew, TAV		2	ST	B02T01
83	04.334.632S	CCHS - Ø 5.5, L 32mm, Langes Gew, TAV		2	ST	B02T01
84	04.334.634S	CCHS - Ø 5.5, L 34mm, Langes Gew, TAV		2	ST	B02T01
85	04.334.636S	CCHS - Ø 5.5, L 36mm, Langes Gew, TAV		2	ST	B02T01
86	04.334.638S	CCHS - Ø 5.5, L 38mm, Langes Gew, TAV		2	ST	B02T01
87	04.334.640S	CCHS - Ø 5.5, L 40mm, Langes Gew, TAV		2	ST	B02T01
88	04.334.642S	CCHS - Ø 5.5, L 42mm, Langes Gew, TAV		2	ST	B02T01
89	04.334.644S	CCHS - Ø 5.5, L 44mm, Langes Gew, TAV		2	ST	B02T01
90	04.334.646S	CCHS - Ø 5.5, L 46mm, Langes Gew, TAV		2	ST	B02T01
91	04.334.648S	CCHS - Ø 5.5, L 48mm, Langes Gew, TAV		2	ST	B02T01
92	04.334.650S	CCHS - Ø 5.5, L 50mm, Langes Gew, TAV		2	ST	B02T01
93	04.334.655S	CCHS - Ø 5.5, L 55mm, Langes Gew, TAV		2	ST	B02T01
94	04.334.660S	CCHS - Ø 5.5, L 60mm, Langes Gew, TAV		2	ST	B02T01
95	04.334.665S	CCHS - Ø 5.5, L 65mm, Langes Gew, TAV		2	ST	B02T01
96	04.334.670S	CCHS - Ø 5.5, L 70mm, Langes Gew, TAV		2	ST	B02T01
97	04.334.675S	CCHS - Ø 5.5, L 75mm, Langes Gew, TAV		2	ST	B02T01
98	04.334.680S	CCHS - Ø 5.5, L 80mm, Langes Gew, TAV		2	ST	B02T01
99	04.334.685S	CCHS - Ø 5.5, L 85mm, Langes Gew, TAV		2	ST	B02T01
100	04.334.690S	CCHS - Ø 5.5, L 90mm, Langes Gew, TAV		2	ST	B02T01
101	04.334.695S	CCHS - Ø 5.5, L 95mm, Langes Gew, TAV		2	ST	B02T01
102	04.334.600S	CCHS - Ø 5.5, L 100mm, Langes Gew, TAV	2	ST	B02T01	
103	04.334.601S	CCHS - Ø 5.5, L 105mm, Langes Gew, TAV	2	ST	B02T01	
104	04.334.602S	CCHS - Ø 5.5, L 110mm, Langes Gew, TAV	2	ST	B02T01	

L004344S

Db. kopfl. KomprSchr. Ø4.5 & 5.5 Steril

B01T01



B02T01



End of document