


PUREVUE™ 4K

4K Direct Coupled & Eye - Piece Arthroscopes



4K Direct coupled & Eye-piece Arthroscopes

DePuy Synthes 4K Arthroscopes incorporate precision machined sapphire lenses with an advanced optical coating, which optimizes image coloration, and enables high quality visualization. This 4K optic design increases MTF (Modulation Transfer Function), and Transmission, while reducing visualization challenges of distortion and image runout in the Arthroscopes. The 4K Arthroscopes equip the surgeon with a clear and high resolution view in order to achieve the best surgical outcomes for their patients.¹



Medos International Sarl
Chemin Blanc 38
Le Locle, CH-2400
Switzerland

© DePuy Synthes 2021. All rights reserved

Please refer to instructions for use, for complete list of indications, contraindications, warnings and precautions.

The third party trademarks used herein are the trademarks of their respective owners.

1- DePuy Synthes, Verification test reports, 2021, ADAPTIV - 103378313 Henke-Sass Wolf Mitek Locking,

103795052 Sapphire 4.0 and 3.5 Design Verification Report

 **DePuy Synthes**
THE ORTHOPAEDICS COMPANY OF 